

Product datasheet for **MC214384**

H2bc3 (NM_175664) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: H2bc3 (NM_175664) Mouse Untagged Clone
Tag: Tag Free
Symbol: H2bc3
Synonyms: H2b-143; Hist1h2; Hist1h2bb
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC214384 representing NM_175664
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGCCAGAGCCTTCTAAGTCTGCACCAGCCCCTAAGAAGGGATCTAAGAAAGCCATCTCTAAGGCGCAGA
AGAAGGATGGCAAGAAGCGCAAGCGCAGCCGCAAGGAGAGCTATTCGGTGTACGTGTACAAGGTGCTGAA
GCAAGTGCACCCGACACCGGCATCTCCTCCAAGGCCATGGGCATCATGAACTCGTTCGTGAACGACATC
TTCGAGCGCATCGCGAGCGAGGCGTCCCGCCTGGCGCATTACAACAAGCGCTCGACCATCAGTCCCGGG
AGATCCAGACGGCCGTGCGCCTGCTGCTGCCCGGGGAGCTGGCCAAGCACGCCGTGTCGGAGGGCACTAA
GGCCGTACCAAGTACACCAGTCCAAGTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_175664

Insert Size: 381 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online >](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_175664.3](#), [NP_783595.1](#)

RefSeq Size: 490 bp

RefSeq ORF: 381 bp

Locus ID: 319178

UniProt ID: [Q64475](#)

Cytogenetics: 13 A3.1

Gene Summary: Histones are basic nuclear proteins that are responsible for the nucleosome structure of the chromosomal fiber in eukaryotes. Nucleosomes consist of approximately 146 bp of DNA wrapped around a histone octamer composed of pairs of each of the four core histones (H2A, H2B, H3, and H4). The chromatin fiber is further compacted through the interaction of a linker histone, H1, with the DNA between the nucleosomes to form higher order chromatin structures. This gene is intronless and encodes a replication-dependent histone that is a member of the histone H2B family. Transcripts from this gene lack polyA tails; instead, they contain a palindromic termination element. [provided by RefSeq, Aug 2015]