

## Product datasheet for MC214419

### Fsd1l (NM\_176966) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Fsd1l (NM\_176966) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Fsd1l  
**Synonyms:** A230072O16Rik; Ccdc10; Csdudf1; Fsd1nl; OTTMUSG00000007002  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC214419 representing NM\_176966  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCCGCGATCGCC

ATGGACTCGCAGAAAGAAGCTCTACAAAGAATCATTTCAACTCTGGCAAATAAGAGTGATGAGATTCAG  
 ACTTCATTGATACACTAAATCACACACTGAAAGGAGTTCAGGAAAATTCTCCAACATCCTTCAGAGTT  
 GGATGAAGAATTTGATAGCTTGTACTCTATATTGGATGATGTAAGGAAAGCATGATTAGCACCATCAAG  
 CAGGAGCAAGTGCGCAAGTCCCAGGAGTTACAGAGTCAGCTCAGTCAATGTAACAATGCCCTGGAGAACT  
 CTGAGGAACTACTAGAATTTGCAACAAGGTCATTAGATATTAAGAACCAGGAAAGATTTTCAAAGGTAAA  
 CAAAACAAAAAGTAAACGAACAGACAAATCAAAAACCCAATTC TAG

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_176966  
**Insert Size:** 396 bp  
**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).  
**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_176966.4](#), [NP\\_795940.1](#)

**RefSeq Size:** 2198 bp

**RefSeq ORF:** 396 bp

**Locus ID:** 319636

**UniProt ID:** [Q8BYN5](#)

**Cytogenetics:** 4 B2