

## Product datasheet for MC214449

### Lrrc58 (NM\_177093) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Lrrc58 (NM\_177093) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Lrrc58  
**Synonyms:** 1810012N18Rik; C330018J07Rik  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC214449 representing NM\_177093  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGGAGGAGGCCGCTCGGGGAGGCCAACTCAACTGGTCCCGCTCAGCGTGTGGCCGAGGCGTTGG  
AGTCGGAGCTGGAGGCGCGCCGAGGAGCGGCGCGCGCAGGCTCTCCTGCGGTGCTGTGCC  
CTACAACCGCTGACGTGCTGCCGCGGGCTGGGCGGGCTTCCCGACCTCCAGCTCTGGACGTG  
AGCGGCAACTCGCTCACGGCGCTGGGGCCGAGCTGCTGACGCTGAGCGGCTGCGCACGCTGCTGCCA  
GGAACAACCGGCTCGGGCGCCGGGCTCGCTGCCAAGGGCTGGCCAGTCGCCGCTCTGCCGAGCCT  
CCAAGTGTCAACCTCAGCGGCAACTGCTTCCAGGAGTGCCCGCTCGCTGCTGGAGTGCAGCGCCTT  
CAGACCTCAGCTGGGCGGCAACCAGCTGCAGAGCATCCCCGCGAGATCGAGAATTGCGGAGTTAG  
AATGTTTATACCTTGGAGGAACTTCATTAAGAAATCCCACCAGAATTAGCAAATCTGCCTTCTGAA  
TTACTTGGTGTGTGTGACAACAAGATCCAAAGTGTGCCTCCCAGCTCTCACAGTTGCATTGCTTGA  
TCCCTCAGTCTGCACAATAACTTGCTGACTTACCTGCCTCGAGAGATCCTAATCTTATCACCTGGAAG  
AGCTGAGTTACGAGGAAACCGTTGGTTGTTTGTAGAGATTTAACGTATGACCTCCGACTCT  
TCTGGAATTAGCTGCAAGGACCATTAAGATCCGAAGTATTTCTTACACTCCCTATGACCTTCTGGGAAT  
CTTCTTAGATACTGGGTTCCAGCCAGCAACTGTCCAAACCAAAGTGTGGTGGAGTCTATTTGACTGTT  
GTGTCAGACAAATTAAGTTTGTGGACTTCTGTGGAAAGTACCGTCTCCCCTCATGCACTATCTGTGCTC  
TCCAGAGTGTCTCCCCTGCAGCTCTGCCTCCATAGCTCCACCTCCCAGAGTGAATCTGACTCAGAA  
GATGAAGCCAGTGTGCTGCACACAGAATGCAGAAAGTCTTCTTGGT**TGA**

**ACGGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_177093



<b>Insert Size:</b>	1101 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"> <li>1. Centrifuge at 5,000xg for 5min.</li> <li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li> <li>3. Close the tube and incubate for 10 minutes at room temperature.</li> <li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li> <li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li> </ol>
<b>RefSeq:</b>	<u>NM_177093.3, NP_796067.2</u>
<b>RefSeq Size:</b>	8592 bp
<b>RefSeq ORF:</b>	1101 bp
<b>Locus ID:</b>	320184
<b>UniProt ID:</b>	<u>Q3UGP9</u>
<b>Cytogenetics:</b>	16 B3