

Product datasheet for MC214496

Cdv3 (NM_175833) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Cdv3 (NM_175833) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Cdv3
Synonyms:	2510010F10Rik; C79446; C230084J24Rik; TPP36
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC214496 representing NM_175833 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCCGAGACCGAGGAGCGGAGCCTGGACAACCTCTTCGCCAAGAGGGACAAGAAGAAGAAGGAAC
GGAGCAGTCGCGCGGCCAACCGCGCGAGCGCGCGGGCGGGAGCAGCGCGCGCGGGCTCCAGGCCGG
CGATGGCGGGAGCCTGGGCTCCGGTGCCAGGTCGGGCGACGGCGGGAGCTTGGGCTCCGGCTCCAGGTCC
GGCAGCGCGGGAGCTCGGGCTCTGGTGTAGGTCGGGAGACGGCGGGAGCTCGAGGTCGGAGACGGCC
GGAGCGCGGGCCCGCCGCAAGGCCATAACGAAGGATGAAAATGAGTGGAAAGAATTTGAGCAAAGAGA
GGTTGATTACAGCGGATTAAGAGTTCAGGCAATGCAAATAGTAAAAGGAAGATGATGATAATGAGAAG
AGAGAAGATCCAGGAGATAATTGGGAAGAAGGTGGAGGTGGCAGTGGAGCAGAAAAATCTTCAGGTCCT
GGAATAAAACCGCTCCGGTACAAGCGCCTCCTGCTCCAGTAACAGTTACAGAAACCCCGGAGCCAGCAAT
GCCCAGTGGTGTATACAGGCCTCCTGGGCAAGGCTAACCACAACAAGGAAAACGCCACAAGGACCACCA
GAAATATACAGTGACACGCAGTTCATCCCTGCAGTCCACTGCCAAGCATGTAGAGAGCCGGAAGGATA
AAGAAATGGAGAAGAGCTTTGAAGTAGTAAGACACAAAAATAGAGATAGGGAGGAGGTTTCAAAAAACCA
GGCCCTTAAACTTCAGCTAGACAACCAGTATGCTGTGCTTGAGAAATCAGAAAGTACAGCCACACACAGTAC
AGCTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_175833
Insert Size:	846 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_175833.2</u> , <u>NP_787027.1</u>
RefSeq Size:	3335 bp
RefSeq ORF:	846 bp
Locus ID:	321022
UniProt ID:	<u>Q4VAA2</u>
Cytogenetics:	9 F1