

## Product datasheet for MC214497

### Cdv3 (NM\_001134426) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Cdv3 (NM_001134426) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Cdv3
Synonyms:	2510010F10Rik; C79446; C230084J24Rik; TPP36
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC214497 representing NM_001134426 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGCCGAGACCGAGGAGCGGAGCCTGGACAACCTCTTCGCCAAGAGGGACAAGAAGAAGAAGGAAC  
GGAGCAGTCGCGCGGCCAACCGCGCGAGCGCGCGGGCGGGAGCAGCGCGCGCGGGCTCCAGGCCGG  
CGATGGCGGGAGCCTGGGCTCCGGTGCCAGGTCGGGCGACGGCGGGAGCTTGGGCTCCGGCTCCAGGTCC  
GGCAGCGCGGGAGCTCGGGCTCTGGTGTAGGTCGGGAGACGGCGGGAGCTCGAGGTCGGAGACGGCC  
GGAGCGCGGGCCCGCCGCAAGGCCATAACGAAGGATGAAAATGAGTGGAAAGAATTTAGCAAAGAGA  
GGTTGATTACAGCGGATTAAGAGTTCAGGCAATGCAAATAGTAAAAGGAAGATGATGATAATGAGAAG  
AGAGAAGATCCAGGAGATAATTGGGAAGAAGGTGGAGGTGGCAGTGGAGCAGAAAAATCTTCAGGTCCT  
GGAATAAAACCGCTCCGGTACAAGCGCCTCCTGTCCAGTAACAGTTACAGAAACCCCGAGCCAGCAAT  
GCCCAGTGGTGTATACAGGCCTCCTGGGCAAGGCTAACCACAACAAGGAAAACGCCACAAGGACCACCA  
GAAATATACAGTGACACGCAGTTCATCCCTGCAGTCCACTGCCAAGCATGTAGAGAGCCGGAAGTACT  
TAAATGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001134426
Insert Size:	708 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001134426.1</a></u> , <u><a href="#">NP_001127898.1</a></u>
<b>RefSeq Size:</b>	4032 bp
<b>RefSeq ORF:</b>	708 bp
<b>Locus ID:</b>	321022
<b>Cytogenetics:</b>	9 F1