

## Product datasheet for MC214536

### Prss34 (NM\_178372) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Prss34 (NM\_178372) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Prss34  
**Synonyms:** C130091P12; Mcp11; mMcp-11  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC214536 representing NM\_178372  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGC**C

ATGTGCTTGGGGATGCTCTGGCTCCTGTTCTCAGTCTCCCTGCCTGGGGAACACGATGCCTCTGACTC  
TTGACCTTGGCTCAGGACAGGGTTGGTAGGCATCGTCGGGGCTGCCCTGTCTCAGCCAGCCGGTTCCC  
GTGGCAGGTCAGCTTGGGCTCTACGACATGGAGCACAGCAGATGGGAACACGAGTGTGGCGGTTCCCTC  
ATCCATCCTCAGTGGGTGCTAACTGCTGCTCAGTGTGTGAGGCCGAAAGAGGTGGAGGCCTATGGAGTTA  
GAGTCCAAGTTGGGCAGCTAAGGCTTTATGAAAACGACCAGCTGATGAAAGTGGTCAAGATCATCCGCCA  
CCCCAAGTTTGTGAGAAGCTGTCTGCCGTGGGGTGCAGGACATTGCTCTGCTGAAACTAGACACCCGA  
GTTGTGCTGTCTGAGCATGTCTACCCGTTTCCCTCCCTGCTGCCTCCCTGAGGATCTCCTTAAGAAGA  
CCTGCTGGGTGGCTGGCTGGGGTGTATCGAGAACTACATGCCACTGCCCCACCCTACCACCTGAGGGA  
AGTGGCGGTCCCATCGTGGAAAACAATGACTGTGAGCAGAAGTACCAGACCAATTCTTCTCGGACAGC  
ACTACCAGGATCATCAAAGACGACATGCTTTGTGCTGGGAAGGAGGCAGGGACTTTGCAAGGCTGACT  
CTGGGGGCCATTGGTATGCAGGTGGAAGTCTTCTGGGTCCAGTGGGTGGTGTGGTCAAGTGGGCATAGG  
CTGCGGGCTTCTGACTTCCCTGGTGTATACCCGGGTGATGAGCTACGTGCTCCTGGATCAAATGCTAT  
GTCCCCACGTTCTAGAACCTTTGAAGGGCCAGATGGGATTCATACAACAGAGGAGACCCAGCCACTC  
ATATTGATCCATTCACTCCTCTCTGCTCCAGCCTCACTTTCT**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_178372  
**Insert Size:** 957 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_178372.2</a></u> , <u><a href="#">NP_848459.1</a></u>
<b>RefSeq Size:</b>	1216 bp
<b>RefSeq ORF:</b>	957 bp
<b>Locus ID:</b>	328780
<b>Cytogenetics:</b>	17 12.53 cM