

Product datasheet for **MC214542**

9830107B12Rik (NM_001177897) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: 9830107B12Rik (NM_001177897) Mouse Untagged Clone
Tag: Tag Free
Symbol: 9830107B12Rik
Synonyms: Trem5
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC214542 representing NM_001177897
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**GCGATCGCC**

ATGGCCTGGGAGCCACATACCTGCTCTCCCCAGTGCTGCTGCTGCTCCTGGCCTCAGGTCCGAGAAGC
 CAAGATTCTCCATCCTACAGTCTTCTAAGTTCACTTCTTCACTGTCACCATGAGTAAGCTCAGGATGAG
 TGACTCGGCATCTATCACTGTGGGATCGCTGCAAATACCAGGATAACTTATCTCAGAGGTATCCAAGT
 GTGGTGTCAAAGATGATTCAAGTCTCCTGCCAATGGGACATTTCCAAGTCTGAGAGCTCACACTTTTA
 GGACCACACAGGACCCGATCTCCATTGCTGCCAGGGGCTCTGTATCCAGTGTTACTATCCCGTGGTCTG
 TGGACTCCTCAGCAAGACTGGTATTCACGGTCTTATTCATTGTACACAGAAGTCATTGGACGACAG
 GCCATGAAGGCACACAACAGTAACCTCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001177897

Insert Size: 450 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001177897.1, NP_001171368.1</u>
RefSeq Size:	3583 bp
RefSeq ORF:	450 bp
Locus ID:	328829
Cytogenetics:	17 C