

Product datasheet for MC214552

Fam163a (NM_177838) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Fam163a (NM_177838) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Fam163a
Synonyms:	A230106N23Rik
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC214552 representing NM_177838 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACAGCGGGAACGGTTGTGATCACTGGCGGAATCCTAGCTACAGTGATCCTCCTGTATCATTGCTG
TCCTGTGCTACTGTAGGCTCCAGTATTACTGCTGCAAGAAGGGCACAGATGGCGAGGATGCTGAGGAGGA
AGAGGAAGAGGAGGAACACGGCCTTTCCATCCATCCCCGAGTCCCCGCTGCAATGCCTGCAGCTCCCAC
GTCTGGACGGCAGAGGGCCCTGGCACCTCTACCAGCGAGTCTGCAGCCAGCCGTGTGGGTGGCCA
GCCACTGTACCACTTGCTCCCCTTACCGCACCCCTTTTACATACGGACAGCTGACATGGTGCCCAACGG
GGGTGGAGGCGAGAGGCTCTCCTTCGCCCTACACACTACAAAGAGGGGGGAACCCCATCCCTCAAGTTG
GCAGCACACAGAATTACCCGGTGACCTGGCCAAGTTCTGGGCATGAGGCCTTACTAATCCAAGGGCTA
TTAGTACCGATGT**A**A

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_177838
Insert Size:	507 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_177838.3</u> , <u>NP_808506.1</u>
RefSeq Size:	3892 bp
RefSeq ORF:	507 bp
Locus ID:	329274
UniProt ID:	<u>Q8CAA5</u>
Cytogenetics:	1 G3