

## Product datasheet for **MC214564**

### Tafa3 (NM\_183224) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Tafa3 (NM\_183224) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Tafa3  
**Synonyms:** 7530404M11Rik; Fam19; Fam19a3; Sa; Sam3; Taf; Tafa; Tafa-3  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC214564 representing NM\_183224  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

**ATGGAGAGGCCACCAGCAACTGGAGCGCAGGCAGCTGGGTGCTTGCCTGCCTCGCTGGCTGTGGA**  
**CGTGTCCGGCCTCTGCTTCCTTGCAGCTCCAACATCCGAGTCTAGTGAAGCAGGGCACCTGCGAGGT**  
**GATTGCTGCACATCGCTGCTGCAACCGGAACCGCATTGAGGAGCGCTCCAGACGGTCAAATGCTCCTGC**  
**CTGTCCGGCCAGGTGGCTGGCACCAGAGCAAAGCCCTCCTGCGTGGACGCCTCCATCGTCTGCAGA**  
**AGTGGTGGTGCAGATGGAGCCCTGCCTGCTGGGAGAGGAGTGTAAAGGTGCTCCCGACCTGTCAGGGTG**  
**GAGCTGCAGCAGTGGACACAAGGTCAAACCACCAAGGTACACGGTAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_183224  
**Insert Size:** 399 bp  
**OTI Disclaimer:** Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).  
**Components:** The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



**Reconstitution Method:**

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

**RefSeq:** [NM\\_183224.3](#), [NP\\_899047.1](#)

**RefSeq Size:** 4317 bp

**RefSeq ORF:** 399 bp

**Locus ID:** 329731

**UniProt ID:** [Q7TPG6](#)

**Cytogenetics:** 3 F2.2

**Gene Summary:** This gene is a member of the Tafa family which is composed of five highly homologous genes that encode small secreted proteins. These proteins contain conserved cysteine residues at fixed positions, and are distantly related to MIP-1alpha, a member of the CC-chemokine family. The Tafa proteins are predominantly expressed in specific regions of the brain, and are postulated to function as brain-specific chemokines or neurokines, that act as regulators of immune and nervous cells. [provided by RefSeq, Jul 2008]