

## Product datasheet for MC214571

### Skint2 (NM\_177865) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Skint2 (NM\_177865) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Skint2  
**Synonyms:** A430070M19; B7S3; OTTMUSG00000008540; RP23-170H12.1  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC214571 representing NM\_177865  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGGAGCTACAGGAGTACTTCTCTGTGTGTCTTACACTTCTCCAGATGGTAACACAGAGTTCAGAAA  
AATTCACAGTGACTGGCTTACAGAGGCCAGTCTTGGCACCATTGGGTGGAATGTTGAAGTGAAGTGTCA  
GTTGTCTCGCCACAACAAGCACAGCATATGGAGATTCGTTGGTTTCGAAATCGTATAGAGAGCCTGTG  
TACCTGTACAGGAATGGTAAAGACCTGCATGGAGAACTATCTCCAAGTATGTGGAGCGGACTGAACTCC  
TGAAAGATGACATTGAAAGGGGAAGGTGACCCCTCAGGATCTTTAACTGACAGCTGATGATGATGGGTC  
ATACCACTGTGTCTTCAAAGTTGGAGAGTCTATGAAGAGCATATCACAGAAATCAAGGTCACAGCCACA  
AGCTCAGTCATGTACATTCTAATGCAGCCCCCTAATATCAAAGGTGTGATGTTGGAGTGTCAATCAGGAG  
GTTGGTCCCACAGCCTCATATGGAATGGAGAGACAACAAGGGAATATCATTCCAGCTACATCAAAGGC  
CCATTCACAGGATGAAAACAAATGTTCAACATGACAATGACCCTTCTATTGAAGCCAGCTCCCACAGG  
AGTATCACTTGTATCTTCAAACCTTCTAACTACCAAGAAGAAAGCATAAGTATTGTCTATCAGGTG  
AACTATTTTCATGAAACGTGTTTGGATAATGATTCTGACTACAATAGGATTCATGATAGCCTTCTG  
CATGACTATTGTGTTTCAGCAACATCTCTTATATGGTACCTTCTTAAGGGGAAATGTCATTGGTTAAAG  
AGTACAATGATATTCATGTTTTTCAGTGATTGCAGTTACAGGAGTCATGCTCATTTTGCATCTAAAACAAA  
GAGGTCTTATGAAATTATTTATGAATGCAGAGAA**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** Sgfl-MluI  
**ACCN:** NM\_177865  
**Insert Size:** 948 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_177865.4</a> , <a href="#">NP_808533.1</a>
<b>RefSeq Size:</b>	1678 bp
<b>RefSeq ORF:</b>	948 bp
<b>Locus ID:</b>	329919
<b>Cytogenetics:</b>	4
<b>Gene Summary:</b>	May act by engaging a cell surface molecule on immature T-cells in the embryonic thymus. [UniProtKB/Swiss-Prot Function]