

Product datasheet for **MC214580**

Fam71f1 (NM_207258) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Fam71f1 (NM_207258) Mouse Untagged Clone
Tag: Tag Free
Symbol: Fam71f1
Synonyms: GARI-L1; Gm764; NYD-SP18
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC214580 representing NM_207258
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGC**C

ATGATGACATCAGTTCACCTAGAAAGTCTTGGTGGACATCGAAGAAGACAGTAAAAGTGATAAGATCAT
ATCCAACCTTCCTTCCCTTGAATGCCTGGGAGAAGTTCAGGGGCTCTTGCTGTGGATGGGAAACGAA
CCCTGGAGTGGGCTGGGAGTGGAGGAAGGACTACTCTGCCAGATGGTTCAATTCAGAAATCAATCTG
TTTCCTAACTCCGTGGTGTTCGAAAGCAACTTTGTTTCAGGTGAGAAGGAGCAGAACTGGAAAGAAATCT
ATAAGGCTTCCAAACAATGGCCCTTGGGGTGACCTCCTCCGTACCTTGCTTCCCTTCCGAACATCCT
CCTTATGGCCAGAGTCAAATGGCATCAGGGGCAGAGCCAGACATGGAACAGACATCCAGAGCGCCAAGC
ATCAACCTGAAGAGCATTCTCCCTTTGAAGTTGTGGAGCTCCAGATCTGGGACGACCAGGAACGCGTCC
TGCGGTTGAGGACAGTCAACGAGAAAATCTACTATCTGAAGCTCCATCCCGACCACCTGAGACTGTCTT
CCACTTTTGGATCCGCTGGTCAAATCTGCATAAGGGTCTGTCCATCACCACCAAGGACCTACCATC
CTTGTCACCTCACTGCCTGGTACCAAGAGCCTTTGCAGCCAGGTGAAAGACTGAGTTAGTGCAGAAGA
AATCGAAAGGCCTCCAGCCAGTGAGAGCCTCACACATCTGATGGCCAGGGGAGAGCCAGACCTCTC
CCAGATCTTTTCCGACTTGATCAGCAGAAAACAGTACAGGAGCGAAAAGATGCACATCAACAAGACCAGC
TCAGAGAAAGCTACTCCCTGTGAAGACAGCATCCCTGTACCTGTGATCTCAACTGGAGGATGCATTCA
TGTTCCGAGAATGGGAAAGAGAAAACCTTCTGGACCAGCCCTCTCACTCCTCGGTACTCTGGCTGC
CTCCAGCAGGCCACGGCTGTCCCTAACTGGAGGAACTCTATA**TGA**

AG**GCGACCG**ACGCGTACGCGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCC
TGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-RsrII
ACCN: NM_207258



[View online »](#)

Insert Size:	1026 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_207258.3 , NP_997141.2
RefSeq Size:	1587 bp
RefSeq ORF:	1026 bp
Locus ID:	330277
UniProt ID:	Q3UZD7
Cytogenetics:	6 A3.3