

Product datasheet for MC214618

Tmem211 (NM_001033428) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tmem211 (NM_001033428) Mouse Untagged Clone
Tag: Tag Free
Symbol: Tmem211
Synonyms: Gm854
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC214618 representing NM_001033428
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCCGCGATCGCC

ATGGAGGTCGGTGTGAGGGCTGCTCTGGGACTCTCTCACCGGCATGTCTGGCCTCAGCCTCATCTCCC
 CCGCCTGGTTCAGAGCCTTTCCTCTTATGGTATCTTACCTACTGCTCTGGCCACCAGGTGACTG
 CTGGAACCAAAGCTGTGTGACCTTTGGGTCATTGAAGGACATGCCAGCTTGGCCTGGAAGTCTCAGT
 GCAACTTCTGGGAGGCTGGCTTTTGCTGTCAATCAGCGCTCTTCTCCTCTCAGCTGGGCACTGGCCC
 CCCGAAGTTGTTTGCCAGGACGGGCTTTGGTCCGACACCGGTGGTGCAGCGGCAGCAGAGTCTCCAT
 GCTTGTAGGTCTGCTGGTTTTCCAGCCACTCTGGCTTCCCAATTTGCCAAAGAAATCTGTGAAGGCTCC
 TCCATGTACCACAGCGGAGCATGCCGGCTGAGCTGGGGCTATGCCACTGCTATCCTCAACATGGTCTGA
 CCAGCTGCTTGTGTCATCAGGTGGCAAGTGACCACAGTCCAGAGAGCAGCCGTCCTTCTCCAGTGAGAC
 TCAGAGAGCCATCCTTGTGCCAGAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001033428
Insert Size: 588 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001033428.2, NP_001028600.1</u>
RefSeq Size:	1827 bp
RefSeq ORF:	588 bp
Locus ID:	333048
UniProt ID:	<u>Q3UUA0</u>
Cytogenetics:	5 F