

Product datasheet for **MC214676**

Slirp (NM_026958) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Slirp (NM_026958) Mouse Untagged Clone
Tag: Tag Free
Symbol: Slirp
Synonyms: 1810035L17Rik
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC214676 representing NM_026958
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGCAGCCTCCGCCATAAAGGGGCTAAGCGCACTGCGTTCCAGTACTGGTCGGCCTATTGCGTTTGTCAGAAAAATTCCTTGACCGCGGCAGCGAGTGAGCTGAGAGAACATTTTGCACAGTTTGCCATGTAAGAAGGTGCACGTACCTTTTGACAAAGAGACTGGCTTTCACAGAGGCATGGGTTGGGTTCAAGTTTCTCTCAGGAAGAATTCAGAATGCACTACAACAAGAACATCATATTGATGGAGTAAAAATCCACGTTCAAGCTCAAAGAGCAAAGGCTTTGCACGGAGCTCAAACATCTGATGAAGAAAGATTCTGAGATGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATTACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-MluI
ACCN: NM_026958
Insert Size: 339 bp
OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_026958.3](#), [NP_081234.2](#)

RefSeq Size: 427 bp

RefSeq ORF: 339 bp

Locus ID: 380773

UniProt ID: [Q9D8T7](#)

Cytogenetics: 12 D2

Gene Summary: RNA-binding protein that acts as a nuclear receptor corepressor. Probably acts by binding the SRA RNA, and repressing the SRA-mediated nuclear receptor coactivation. Binds the STR7 loop of SRA RNA. Also able to repress glucocorticoid (GR), androgen (AR), thyroid (TR) and VDR-mediated transactivation (By similarity).[UniProtKB/Swiss-Prot Function]