

Product datasheet for **MC214739**

6430531B16Rik (NM_001033465) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: 6430531B16Rik (NM_001033465) Mouse Untagged Clone
Tag: Tag Free
Symbol: 6430531B16Rik
Synonyms: Gm1093
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC214739 representing NM_001033465
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGCTCAGTGCTCCAAGCCGAGCTCAGCTTCGTTCTGCGCAACGCCTGGGGATATTCATAACCTCTGCG
CCTGGCTGGACCGCTCCCACTCAGCCGCCCAAGCGCCACTTGGCTCGGGACTTCAGTGATGGTGTGCT
GGTGGCCGAGATTGTGAAGCACTCCGCCCTCGGCTGGTGGACTTGCACAGCTACGTCCCAGCCTGCAGT
ACAGACCAGAAGCTGAGCAACTGGAGCCTTCTCAACAGGAAAGTCTTTCGCAAGCTGCAGTGTGTATCT
CAGATAATGATATCCAAAGGGTGGTATCCAACAGGCCTGGAGTTATCGAGTCTATTCTATGTGCACTGAG
GGAGAAGATGGAGGCCAGGACTGTTATGTAGGCTCAGCTGGTACAGCTGATCCAAGACACTCTGCTGTG
GATGCTGCCGACCTGGCGGGATCTTACCCCCAACCCACATACCTCTAAGATGCAAGATCCTGATGTGG
ATAACAGCAGCTTGTGAGAACCCTCTGGGACTCCGGTGGAGAGCACCCTCTGTTCCCGAAAGCTGGGT
CCATGATGTTGAACCCCAAACCTTCAAAGAAGGCAGGCCTACAGAGCCCTCCAGCCCCCTCCAGCAGC
AAGATCCTTCAAGATCAAAGGGATTTGGAGAAAATGGGCTGCTGTGCTTGCAGAGGGGACTCAATTGAAG
GGCTCTGGGATCACCTGGTCAGCATCCAGCAGCAGCTGGAGGATAAGGAGCAAGCACTAGCCATTCTGCA
GGAGACGGTCAAGATTTTCAGTTGAAGGTGATGAGGTTGGAACACCTGGTACAGCTGAAGGACCCAGCGG
ATCTGGGAGCTGATGAGAACTGGACCCGAGGAGCATCAGAGCTGGAGGAGGCCCCACCCAGAGTCACTCA
GACCTCGATCT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001033465
Insert Size: 924 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001033465.2</u> , <u>NP_001028637.1</u>
RefSeq Size:	1161 bp
RefSeq ORF:	924 bp
Locus ID:	381933
Cytogenetics:	7 F4