

Product datasheet for MC214745

Pabpn1l (NM_001007462) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Pabpn1l (NM_001007462) Mouse Untagged Clone
Tag: Tag Free
Symbol: Pabpn1l
Synonyms: ePABP2; Gm1108; Pabpn1
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC214745 representing NM_001007462
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGAGCCTTACCTGAGCAATGAGCTCTCCACCCCCCACTGAAGCCTGGCTCCAGACGGTCTCCTCAG
ATCCTGAGGCCAAGGTTGGGGGCTGGGGCAGGACTGAAAAGACCTCCTTGGTGCCAAGGGCTGGGAG
TAGAGCAGGAAGTGACAAGGAAGCAGAGGAAAATGAAGATGCGAGCTTCTTCTCCTCTGCTGGAACCA
GAGAACCCTGGCCAAGTCCCCAGTATTCAACCAGGAGTTGGAGGCCATCAAGCTGAAACTGTGGGCCATGG
AGCATGCTGAAGCCCAGCCGGAGCCACCCTGTGTACAGAGAAAGCCACAGAGGAAGAGAGGGCTGAGGT
CAGACAGCTGCTGAGCCCTGAGACCGTAGACTGCTTCTTCTCTAGGACCTCCAAGGAAAATGTGGAGGCT
GACCACAGATCTGTCTTTGTAGGCAATGTGGATTATGGTGGCTCAGCTGCAGAGCTGGAGGCCTACTTCA
GTCCCTGTGGGGAGATTCACCGGCTCACTATCCTGTGCGACAAGTTCTCTGGACACCCCAAGGGCTATGC
ATATATAGAGTTTGCTTCGCACAGATCTGTCAAGGCTGCTGTGGGGCTGGATGAGAGCACCTCCGAGGT
CGGGTCATCAAGTACTTCCAAAAGGACGAACCTCCAGGCATCAGCTCCACAGACCCGTGGAGGGCTGC
GGACACACTCGGGCAACAGGGCTGCCTTCTGCACGGTAGCCTCCATCGAAAAGCTCGGTTAAGAGCCCA
TGGACGGAGCCGGGCCATGGAGGAGCACACAGTGGTTCTCTCCCTAC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001007462
Insert Size: 822 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001007462.1</u> , <u>NP_001007463.1</u>
RefSeq Size:	1144 bp
RefSeq ORF:	822 bp
Locus ID:	382035
UniProt ID:	<u>Q5XFR0</u>
Cytogenetics:	8 E1
Gene Summary:	Binds the poly(A) tail of mRNA.[UniProtKB/Swiss-Prot Function]