

Product datasheet for MC214763

Rhox12 (NM_001025083) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Rhox12 (NM_001025083) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Rhox12
Synonyms:	Gm1148
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC214763 representing NM_001025083 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGATCGCC**

ATGGCTCTCGAACCCCATCATGTGGACCCCAATTTCTACAACTGGGGGAGAATGAGATCGAAGTGACCC
 TTGATGCTGATCAAGAGGCTGATGGAGCAGCAGAGGGAGGCAGCTTTGGAGACGGTTCTCTAAATGGCTC
 AGACAAATTAAGTGCCAGGGCATCCCAGACGACAAGGATGATGTGATCTACGTTGGAGAATTGAAGAAC
 ATTGGCAATGACATCAAAGACGAGTGCCATGGGAGCCATCAAGGGTCTGGAGATCCACAGCCGGAGGAGA
 AGCAGAAGAACTCGGCGGCTGCCAGAGTCCACAGGTCCGGCGCACAAGGCCACGGATCCAGTTGGGTTT
 CACGCCCAGGCAACTGAATGAACTGGAAGACTTTTTTGAAAAAACTAAGTACCCTGATGCACTCACAAAGG
 AAGAACCCTTGCAAAGCACTTGTACCTAGCAGAATCCAAAGTTCAGAGATGGTTTAAGAAAAGAAGAGCCC
 ATTATAGGAAAGAAGCAGCAGTCTCAAGTGCTGCAGTGTGCATCTGCTGATGGCCAGGATGCCATGCAGTG
 A

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001025083
Insert Size:	561 bp


[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001025083.2</u> , <u>NP_001020254.1</u>
RefSeq Size:	935 bp
RefSeq ORF:	561 bp
Locus ID:	382282
Cytogenetics:	X A3.3