

Product datasheet for **MC214764**

Sly (NM_201530) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Sly (NM_201530) Mouse Untagged Clone
Tag: Tag Free
Symbol: Sly
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC214764 representing NM_201530
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGAAGAAATGGCTCTTAAGAAATTGAAGGTGATACCAAAGGAAGGTTACTTATTACTTTGGACTTTG
ATGATGAGGACGATGACATAAAAGTTTCAGAGGAGGCTCTTTCGGAAGTAAAGAGCCCAGCATTTGATAA
AAATGAGAATATATCGCCTCAAGCAGAAGCAGATGAAGATATGGGAGATGAAGTAGACAGTATGTTGGAT
AAATCTGAAGTAAATAATCCAGCAATTGGTAAAGATGAAAATATATCGCCTCAAGTAAAAGGAGATGAAG
ACATGGGACATGAAGTAGGCAGTATGTTGGATAAATCTGGAGATGACATTATAAGACGCTTCACATAAA
GAGAAAATGGATGAAACTTATGTCAAAGAATCTTCAAAGGCAGCAACCAGAAAATTAGAAAGATTTTGC
AAAACGAACGAACGAGAGAGGAAGAACATCAACAACAAATTTTGTGAGCAGTATATAACTACATTTCAGA
AGTCTGATATGGATGTACAGAAATTCATGAAGAAAAAGAAAAATCAGTGAATAGTTGTCAAAAAGAAC
ACAAGCATTGAACTGTCCAAATGTAGTCAGAACCAGACCCTGGAAGCAGTTAAAGAAATGCATGAGAAG
TCCATGGAGGTTTTGATGAACTTGGGGACCAAGAACTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_201530
Insert Size: 669 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_201530.2</u> , <u>NP_963288.2</u>
RefSeq Size:	918 bp
RefSeq ORF:	669 bp
Locus ID:	382301
Cytogenetics:	Y C2