

## Product datasheet for **MC214807**

### Tex28 (NM\_001126488) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Tex28 (NM\_001126488) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Tex28  
**Synonyms:** Gm1538; Gm1811  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC214807 representing NM\_001126488  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**GCGATCGCC**

ATGGTTTTAAAGGTGAAAGTACCAAGAGCTCAAGTGCAACCCTCCCCACCAACCTGCCCTTTACAGGT  
CTCTGTCTTCTCTGTGAGGATTGCCCTAGTAGTCACACTAGCTTCTCAGATGGTGAGCTTGCCCGCAA  
TGTGAGGGAAGGTGTCAAGCATCGAATCTTCTACCTCTCAGAGCAGTTGAGAGTGGAGAAGGCCAGTAGG  
GATGAAAATACCATGAGCTACCTCAAAGCTGGTATCCAAAGCTGACCGGCACCAGGCCCCACACATCCGGA  
AGGCCCTTTGAGAGGGTGAACCAGCGCACCTCTGCCACTATTGCTCACATTGAACGAAAACCTATCAGTG  
TCATCAGCAGCTGAAGGAGCTGGAAGAGGGCTGCAGTCCCACAAGCTCAGTGCTGAAAGTGGGCAGTGGC  
TTGGACAGCCATAAGCAGCCAAGTGGGAAAGTTTCATATTCCAAATTGTCCAAGCCAGGTGGGGAAGATA  
GCCTGCCCATCAACGTTGCTAGGTCCTCTACTCTAGAGAGTCACTTATCAGAAATGCAACAGAGGAAATT  
CTCAGACAAAAGTATGTGGCCAGCAACAAAAGCTGCTTTTGCAGAAGATGAAGGAAGAACTGACAGAA  
GCCAAGAAGGTCCATGCTAGCTTTCAGGTATCCCATCAAAGCCTCAAAGAAAGTCATATGATTGATGTGC  
AGAGGATACTGGAGTCACTCCAGGAAAAGAAAACCAACAATCACTGATGGAAAAACAGGTGAATGATCA  
CTTGACAGAGATACCTGGATGAGATTTGCCACCTCAAACAGCATCTGGCATGCACTGAAGAGAAAATGGCC  
TATTTGTCCTATGAAAGAGCCAAGGAAATATGGGATGTAATGGAGATTTTCAAGAGCAGAATAACCAAGC  
TAGAAACTCTACAGCAAGCTACCAACTGGAGATGATGGCCAGCCTCAGAACCCGTCCTCAAGGACTTTTT  
GTTGAGATTCATAAGCCTGCTTCTCACTCTGACCACTATCCTCTTGGTCGTGGTCTCTACATTGTGTTCT  
TGCCCTTGCCCTTGCTCAGCTCACGCCTGCGCATCTTATTGATTTATGATCATTGGGCTTGGGACTT  
TGGCCTGGCAGAAGCGCCATGTCATCTCTATTATTGACTGGCAGGCTGGTCCCCTTTAAGTGGAGACA  
AGACTTAAAGGATGCTAAGCCTCCATCAGATGGACAC**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

<b>Restriction Sites:</b>	Sgfl-Mlul
<b>ACCN:</b>	NM_001126488
<b>Insert Size:</b>	1230 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001126488.2</a></u> , <u><a href="#">NP_001119960.1</a></u>
<b>RefSeq Size:</b>	1420 bp
<b>RefSeq ORF:</b>	1230 bp
<b>Locus ID:</b>	385380
<b>Cytogenetics:</b>	X A7.3