

Product datasheet for **MC214849**

H2-M10.6 (NM_201611) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: H2-M10.6 (NM_201611) Mouse Untagged Clone
Tag: Tag Free
Symbol: H2-M10.6
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC214849 representing NM_201611
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGAGGAAGCCTGGACCCTGGACCTCCTCCTCTGTTGATGGCAATAGACCTGACCCAGTATTGTGCAG
GCTCACACTGGCTGCAGATTTTCAACATTGTGATTTTGGAGCCTGGCATGTTGGAGGCCCGTTTCATCCA
GGTCGGATACGTAGACTCCATACAGTTTCAGGGATTTCGACAGCAAAGAACCAAATGCAAGGATGCAACCT
CGGGCTGCTTGGATGAAGCAGGACCCACCTGAATATTGGGAGAAAAGACAGCAAAGCCCTGAGCTTCT
CACTGTACAGGAAGACGAATTCTTAAGTACATGATGTCTGCAAATAAAGACAAGAATAATGGTTATCACAC
CCTGCAGGAAGTGTATGGCTGCAACGTGGCAAATGATGGGAGCTTCTCGAAGGCTACTACCGGCTCACA
TACTATGGATATGATTCTTAATCCTGAATGAGGACCTGAACTCTGGACTGCAGAGGGTGAAGTAGGTG
AAAAATTCAAGACATACCAGGAGCAAGGAGGCATGACAGAAAGTTGGAGGACTTACTGCTAGGGGAATG
CGTGAAAGGCTTCTTAGATGCCTGGACCTTGAAAGAAGACCTGCTGCACTCTGATGCCCCCGAACA
CATGTGACCCACCATGTAAGACCTGAAGGAATGTACCCTGAGGTGCTGGGCCCTGGGTTTCTACCCAG
CTGACATCACCATGACCTGGAAGAGGGATGAAAAAACACACCCAGGACATGGAGCTGCCAGACACCAG
GCCTGCAGGGGATGGAACCTTCCAGAAGTGGGAGCTGTGGTGGTTCTCTGGAGAAGAGCTGAGATAC
ACATGTCATGTCAACCATGAAGGTTACCTGGGCCACTCACTCTGAAATGGGAGCCTCCACACACCATCC
CCATCATAGCAATTCTCATCGGCCTGGTCTTGGAGCTTTGGTGGTGGAACTGTGGTATTTTTCTGGT
ATGGAAGAATTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_201611
Insert Size: 993 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_201611.2</u> , <u>NP_963905.1</u>
RefSeq Size:	1405 bp
RefSeq ORF:	993 bp
Locus ID:	399549
Cytogenetics:	17 B1