

Product datasheet for MC214855

Tdpz5 (NM_207273) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids

Product Name: Tdpz5 (NM_207273) Mouse Untagged Clone

Tag: Tag Free

Symbol: Tdpz5

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >MC214855 representing NM_207273
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGTCAGAGGACATGGAAGTCACGAACTGGGGCTACACACATATCAGTGTTAAGGAATTCTGCTATGTG
GGACTATTAGGAACCTTTACCTTGCATAGATGGAATTAGGAGAACGATTACGAGCCAGTCTTCTCGTT
AGAGGCCAATGATGAAGTGACATGGTGTGAGAGCACACCCAAATGGGGTCGATGAAGTAAGCGAATGT
TACATGTCAGTTTTCTGGAGTTGCTCAGCTGTCGAAAGAGCCAGTTTGGGCAAAGTACGAGTTCTGGA
TCACAACCTCCCAAGGAGAGAAATATCAATGTATGAAGAGCTTCAATGTCCATAGTTTTCAAAAAACCA
ATACAGGGGATTCAAAAAGTTTCCTCTGGAGATTTCTCATCTCTCATCCACGTCGGTTTTCTCCCTGAA
AACAAGCTCACCTCTGCTGCAAGGTGAGCATAGTAGGATCTGTCTTTGGCATGCCTGGACAGAACATAA
CACCTGCAATCAAGGATCCAAGGCACCTGTTGACAGATGACCTAGGAGAGCTGTGGGAGAATTCCTGTT
CACAGACTGCTGCCTGTTGGTAGCTGGCCATGAATTCAGGGCTCACAAGGCCATCCTAGCAGCTCGCTCT
CCAGTTTTAGAGCCATGTTGAACATGAAATGGAGGAGAGACTAGGAAACCCACTGAGATCCATGACC
TGGATCCCAAAGTCTTTAAGGAGATGATGGGCTTCATTTACACTGGGAAGGCGCCACACCTCCAAAGCCA
CTCCATGGCCACTGATGTGCTGACAGCTGCTGACAAGTATGGTCTGGAGGGCCTGAAGTCTGTGTGAG
GATGCCCTCTGCAGGAACCTCTCTGTGGAGAATGCTGCACAACTCTCATCTGGCTGACCTCCACAAA
GAGAGCAGCTGAAGACTCAGGCCCTGTATTTTCATTGCACTTCATGCTTCTGTGGTCTCTGAGACCTCAGA
GTGGAAGTCAATGATGGAGACCCATCCCCACTTGTTGGCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_207273



[View online »](#)

Insert Size:	1023 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_207273.2, NP_997156.2</u>
RefSeq Size:	1023 bp
RefSeq ORF:	1023 bp
Locus ID:	399676
Cytogenetics:	3 F2.1