

## Product datasheet for **MC214859**

### Plcxd1 (NM\_207279) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Plcxd1 (NM\_207279) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Plcxd1  
**Synonyms:** A330045H12Rik; Gm576  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC214859 representing NM\_207279  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGAGCCAGTCAGAGTTGGTCCCCACCCACCCCTAATGCAGGGACTAGAGCCGAGCAGGGTAGGAAGAG  
TGGGTGAAGTCTGCCGCTCAGTGCATAAGAACGTGTTGACCGTCCAGAGCACACAGTAGGCACCGAAGC  
TGCAGCTCCCATGGGAGACCAGGCAGACTGGATGTCGACGCTGTGCCCTCAGCTGTGGGATGTCCCCTG  
CACCACCTGTCCATCCCAGGGAGCCACGACACGATGACATACTGCCTGAACCGTAAGTCTCGGATCTCCC  
GTGCAAGTTCGTGGCTACTGCACCTGTTGGGCCGGTCTGTCATTTCATACCCGGCCCCGTGGTGTGAA  
ATGGTCGGTACACAGACCCTGGACGTGACGCAGCAGCTGGATGCGGGAGTGCAGTACCTGGACCTGCGG  
ATCGCGCATGCGCCGGAGGGCTCTACACGGAACCTATGCTTCGTGCACATGATGTACACGAAGCGCTGG  
TGGAGGACACCCTCACTGAGATAGCTGAGTGGCTGCAATCACACCCCGCGAGGTGGTCACTTGGCTTG  
CAGGAACCTCGAGGGGATGACTTGTGAGTTGCACGACTATCTCGCGGGCTGCATCGTGAACATTTTGGG  
GACATGCTGTGTCCAGTGGGAAAGTCCCACACTGAGGCAGCTATGGGCAGGAGCAGCAGGTCATCG  
TGTCTACGAGGATGAGGCCACAGTCAGCCGTTATGATCAGCTGTGGCCTGCGATCCCATTGGTGGGG  
AAATGCAGTGAAGACAGATGTGCTGCTGCGGTTCTGGAGACCATGAAGGGCCAAGGCCGCCAGATGGT  
CTGTTCTGTTGCTGGGATCAACATCACTGAGAACCCTGTGCTACATCCTCCTGCACCCAGTGGACTCCCTGG  
AGGAGATGACACGTCGCAGCCTCCCCTCATGACCAGTGGGTGTGTGCACAGCAGCCGGGGCAGAGCCC  
TCAATGCACCAACATCATCGCAGGCGACTTCGTGGATGCGGATGGCTTCGTGAGTAAGGTCATCAGCCTG  
AATTGCAAACCTGCTGTCTCC**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAAGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_207279



|                               |   |
|-------------------------------|---|
| <b>Insert Size:</b>           | 1074 bp   |
| <b>OTI Disclaimer:</b>        | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).  |
| <b>Components:</b>            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).  |
| <b>Reconstitution Method:</b> | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| <b>RefSeq:</b>                | <u>NM_207279.2, NP_997162.2</u>   |
| <b>RefSeq Size:</b>           | 3484 bp   |
| <b>RefSeq ORF:</b>            | 1074 bp   |
| <b>Locus ID:</b>              | 403178  |
| <b>Cytogenetics:</b>          | 5 F   |