

Product datasheet for MC214867

Mrgpra7 (NM_207535) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Mrgpra7 (NM_207535) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Mrgpra7
Synonyms:	MrgA7
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC214867 representing NM_207535 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGACGAAACCAGCCCTAGAAGTATTGACATCGAGTCACTGATCCCAAACCTTGATGATCATCATCTTTG
GACTGGTTGGGCTGACAGGAAATGCCATTGTGCTCTGGCTCCTGGGCTTCTGCTTGACAGGAATGCCTT
CTTAGTCTACATCCTAACTTGGCCCTGGCTGACTTCCTTCTCTCTGTCACTTCATAAATTCAGCA
ATGTTTCTTCTCAAGGTTCTATACCCAACGGTATCTTTGTCTATTGCTTTTACACCATCAAATGGTTC
TCTACATCACAGGCCTGAGCATGCTCAGTGCCATCAGCACTGAGCGCTGCCTTTCTGTCTGTGCCCAT
CTGGTATCACTGCCGCCGCCAGAACACACATCAACTGTCATGTGTGCTGTGATTGGATCTTTTCCGTG
TTGATCTGCATTCTGAAAGAATATTTCTGTGATTCTTTGGTACCAAATTGGGAAATTACTATGTGTGTC
AGGCATCCAACCTCTTTATGGGAGCATACCTAATGTTTTGTTTGTAGTCCTCTGTCTGCCACCCCTGGC
TCTGCTGGCCAGGTTGTTCTGTGGTGCTGAGAAGATGAAATTTACCAGATTATTCGTGACCATCATGCTG
ACCATTTTGGTTTTTCTCTCTGTGGGTTGCCATGGGGCTTCTTCTGGTTCCTGTTAATCTGGATTAAGG
GTGTTTTAGTGTACTAGATTATAGACTTTATTTGGCATCAATTGTCCTAACTGTTGTTAACAGCTGTGC
CAACCCCATCATTTACTTCTTCGTGGGATCATTCAGGCATCGTTGAAGCACCAGACCCTCAAAATGGTT
CTCCAGAGTGCACTGCAGGACACTCCTGAGACACATGAAAACATGGTGGAGATGTCAAGAATCAAAGCAG
AGCAG**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_207535



[View online »](#)

Insert Size:	918 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_207535.1</u> , <u>NP_997418.1</u>
RefSeq Size:	918 bp
RefSeq ORF:	918 bp
Locus ID:	404236
UniProt ID:	<u>Q91ZC5</u>
Gene Summary:	Orphan receptor. May be a receptor for RFamide-family neuropeptides such as NPFF and NPAF, which are analgesic in vivo. May regulate nociceptor function and/or development, including the sensation or modulation of pain (By similarity).[UniProtKB/Swiss-Prot Function]