

Product datasheet for **MC214875**

Vmn1r179 (NM_207545) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Vmn1r179 (NM_207545) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Vmn1r179
Synonyms:	V1rd17
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>MC214875 representing NM_207545 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGC**C

ATGTCCTATCTAATCATCACTCCAGTGTCTGCGAAAGCCCTCTCTATGTGTACGCTTGATAAGTCGCTGA
AAACCACTGAGGAAGTAGCTCTTCAGATCTTCTGCTTTGCCAGTTTGGTGTGGAACTGTGGCCAAATGT
CTTTCTGTTTGTCCATAATTTCTATCCAGTCTTGACTGGTTCTCAACAGAGTCCCAGACAGGTGATTTTA
AGCCACATGGCTGTGGCCAATATCTTGATTCTCTTCATCACTATATTTCCAACCAATGTGATGACTTTTG
CTCCAAGAAATCCTCCAAGTACCTCAAAATGTAATAGAAATCTTCATTCGCCTAGTGGCAAGAAGCAC
AACTTGTGTTCCACCTGTGTCCTGAGTATTTATCAGTTTGTCACTCTTGTTCCTGTTAATAGAGGTAAA
CTCATACTCAGAGCAAGTGTCCCAAGTTTGGCAAGGATTCCTGTTACAGTTGTTGGTCTTCAGTGTCT
TAAACAACATCTACATTCCAATTAAGGTCAGTGGTCCACAGATAGCAGACAATAACACTGACTGTAGAAG
CAAGTTGTTCTGCTTCACTTTTGGATACAGTAGAGTCAATTTCTTCTTGAGGATTGCCTATGATGCCACA
TTCATGAGCATCATGGTCTGGACCAGTGTCTCCATGGTACTTCTCCTCCATAGACATCACCATAAAATGC
AGCATATCTTTACTACCAATCGACTACAGAGGCCAAGCTGAGACTAGAGCAACCATATCCTGAT
GCTGGTGGTACATTTGTTAGCTTTTATCTTCTAAATTTACTTGCATCATCATTCAAGCCCTTGTTACA
GATTCTCATGATTTGTGAGGCATTGGAGTGAAGCTTTGGTACAGGCTTTCCCACTATCTCTCCCTTAC
TGTTGATCTTCAGAAAATCTAAGGATCCTGTTCTGTGTTCTTCCGATTGTTGAAATCCT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_207545
Insert Size:	972 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_207545.1</u> , <u>NP_997428.1</u>
RefSeq Size:	972 bp
RefSeq ORF:	972 bp
Locus ID:	404286
Cytogenetics:	7 A3