

## Product datasheet for **MC214886**

### Olfr250 (NM\_207552) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Olfr250 (NM\_207552) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Olfr250  
**Synonyms:** MOR170-14  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC214886 representing NM\_207552  
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**GCGATCGCC**

ATGCTTCTACAGATGCAACACATGAAACAAATGATTATGGAAAATGACTCTTCAGTGTCTCAGAGTTCATAC  
 TTATGGGACTGACATACCAGCCTGAGCTCTGGTGGCCCTATTTGTTCTGTTCTTGGTGAACATACAGC  
 CACTGTGATGGGAAACCTGAGCTTAATGACTCTCATCTGTCTGAATTCTCATCTGCACACACCAATGTAC  
 TTTTTTATTCTCAACTTGTCTTCAATTGACTTCTGTTATTCATTTGTTTTACCCAAAAATGCTGATGG  
 GATTTGTCTCAGAACACACACCATCTCCTTCACAGGATGCATGACTCAGCTATTTTTTTCTGTCTTTT  
 TGTTAACTCTGAGTGCTATGTGCTGACAGCCATGGCCTATGATAGATATGTGGCCATCTGTAGGCCCTCTG  
 CTGTACACAGTAGTTATGTCTCCAGGGCTTGTCCCTGTTAATGCTTGTGCACACTTGATGGGTGTCT  
 CTTCTGCTGTTGTACACACAGGATGTATAATCAGCTCAGGTTTTGTGGTTCAAAGTAATCAACCACTA  
 CATGTGTGATACTTTCCCCTCCTTGAGCTCTCCTGTGGTAGCAGTCATGTCAACGAGCTTGAAGTTCT  
 GTTTCTGTGGCTGTCGTTGTTGTTATCTAGCCTAATTATTGTATCCTCATATGCTTTGATTCTTGTC  
 ATGTTATCCATTTGTATCATCTAAGGGTTGGTCAAAGCTGTGAGCACATGTAGCTCTCATATAATAAC  
 TGTTGCCCTGTTCTATGGATTTGGTTTGTCTGCTCATATCAAGCCATCATCTGCAGAATCTGTAGTTGAG  
 AGGAAATTTTTTTCAGTAGTTTATACTTTTGTGCTGCCTTTGTTGAATCCGCTCATTACAGCCTCAGGA  
 ACAAGGATGTCAAACCTGCTTTGAAGAGAACACTAAAGACAGTTACAATCCAAGGGAAGTCTTATGTTG  
 TAGCCATAAATCT**TGA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_207552



[View online »](#)

<b>Insert Size:</b>	996 bp
<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<a href="#">NM_207552.2</a> , <a href="#">NP_997435.2</a>
<b>RefSeq Size:</b>	996 bp
<b>RefSeq ORF:</b>	996 bp
<b>Locus ID:</b>	404312
<b>Cytogenetics:</b>	9 A5.1
<b>Gene Summary:</b>	Olfactory receptors interact with odorant molecules in the nose, to initiate a neuronal response that triggers the perception of a smell. The olfactory receptor proteins are members of a large family of G-protein-coupled receptors (GPCR) arising from single coding-exon genes. Olfactory receptors share a 7-transmembrane domain structure with many neurotransmitter and hormone receptors and are responsible for the recognition and G protein-mediated transduction of odorant signals. The olfactory receptor gene family is the largest in the genome. The nomenclature assigned to the olfactory receptor genes and proteins for this organism is independent of other organisms. [provided by RefSeq, Jul 2008]