

Product datasheet for MC214919

BC051142 (NM_001001177) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	BC051142 (NM_001001177) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	BC051142
Synonyms:	NG8; Tesb; Tesbp; TSBP
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC214919 representing NM_001001177 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGGCAGATCTCGAAATAGTTTTGGCTGTATCCTGACTGTATTGGCCCTTGCTATTCTGACTATTCTGA
TAGCAAGATGGTCACGTCGTAGGAACAATGGGGCAGGGATGTCCAGATTCGCTTCAGAACAGAGTATTTCC
CATGGGAGACACAAAAGACGAGGTAAACAGTTTCGTTTCCTATGCAGGACTGGCAGTGTATCAATGAG
GACCTCACTGTGTCTCTTTCTTTAGGATCTCGTTGCTCGTTTCATCAGATAGCGAACTGACTCAAA
AAACCTTCCTGCAAACTCCTGGACCTATTGTACAATGTATGGATTCCAGTGAAAGCTGCACCCCTGAGGG
ATCCAACGCTGTAATACCTGTGGACATATGTTGTCAACCAGCACTGAAGCCTAGCTGTTGTCAACAGCTA
GCAATCAAACCTAGCTGTTGTCAAGCCTTGTCAAGCCTTGTCAAGCCTTATCAACCTTGTCAAGC
CTTCTCCGCCAGACACCTTGCCCTTATTAATTACACAGAGAACTCAACAAAATGTGGACCACCTATAAT
ATTGAAAAGCCATCAGATTGCAAAGGCCAGATGGTAAGACATGGAGGGAGTGAACCTTCCAGGGTATC
GGAGATGACCAATAAAGTGCTGTAAGGAGTGGAGAAGATGCCAGTGTGAATACTGCTGTCCAAAGAGAT
TCGAGCCCCACAGACGAGGAACCATATGGTGTGAACAACCAAGCTCCTTCTCTTCATCAGCATCAGCA
TCAGC
GGTGACCAAAGAATGAAAAGTGAATTTGATGCCAGACCATATGTTATGAAGAGCCATTCCTTTTCAGAA
GACTGGAGTCCCTGATAAGGAGGAACCAAGACGGGATAGTAATAAAGCAGAATGAAGATGGCTTTCTGAT
AAAGAAGGCAGAGGAAGGCAAAAAGAAAAAGAAACATCCAGGAGGGCATCAGCACCCAGAAACAGAT
TCTGATTCCAGCGCTCCACAGACCGGAGAGGCTCTTTTCGAAAGTAAACACCCGAAGAACAACAGAGG
GAGAAGAGAAGAAACGCAAGAAAGAAAGAAAGAAAGAAAGAAAGAAAGAAAGAAAGAAAGAAAGAAAG
GCAAGACGGCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCTCTGGATT
ACAAGGATGACGACGATAAGGTTTAA



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001001177
Insert Size:	1182 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001001177.2, NP_001001177.2</u>
RefSeq Size:	1900 bp
RefSeq ORF:	1182 bp
Locus ID:	407788
Cytogenetics:	17 B1