

Product datasheet for MC214937

Ccdc92b (NM_001013784) Mouse Untagged Clone

Product data:

| | |
|---------------------------|------------------------------------------------------------------------------------|
| Product Type: | Expression Plasmids |
| Product Name: | Ccdc92b (NM_001013784) Mouse Untagged Clone |
| Tag: | Tag Free |
| Symbol: | Ccdc92b |
| Synonyms: | E130309D14Rik |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Fully Sequenced ORF: | >MC214937 representing NM_001013784 Red=Cloning site Blue=ORF Orange=Stop codon |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAGCGCAGCTAGCGCGGGAGAAGAGCTCTTCCAGAGACCCAGCGAGGCCCGGGCCGGCCCGGGCT
CTCTCGGGTCCGAGGTGGGGCCAGGCCCTGGTTGCCATGGATACCGTGTCCCTGGAGCATCAGATCCA
GAGCGTCAACGCCACATCAGCTTCTGAAAAAGGAGCAGATGGCCCTGCTTCGAGACCTGCACCTGGAG
ATCCTGAGGCTGCAGAAACGCTGCTCAGAGCTGACCCACGACCTGGAGATGAGAGAGGTCCAGTCTACC
AACAGAAGAGGGCTCCCGAGAGCTGGAGGACAAGTGCCGCGCGTAGAGTCGCAGCTGGAGGCACGGGC
GGCGGCCAACGAGGAGCTGCGGCCGGGAGGTGGCACAGCGGGAGGCGCTGGTGTCCGCGCTGCGCTGACG
CTGCGCTCGGAGGAGCGCCGCTTCTGGAGGAGCTGCGGCCGCGCAGCCATCGCGCCACGGTGTGGGCA
CCGAGCTGCAGAAGCACACCGAAGCGGCCGCTACCTCTCCTGCCAGCTGCACGCCGCGCGCCAGAGACT
GCAGGGCCCCCGCTCTGGCGCCGCGGCCGCGGAGCCCCGGTCCGCCGCGTGCACAACGGGCCCGCCG
CCTCCGGAGGGCTGCCAAGGGTCCGAGCCGGACTGGGCGACCTGGGAAGACGCCGATCCCATGCCGG
ACCCCGCGCTCTTCTCTACCCACGAGGCCCGCGGCCCGGCCCGTGGTCTCGCCCTCTGCTACG
CCAGGAGATGCCGACCAAGATGCCCGCACCCCGTGCCACACAGGAGCCCTGTACCAAGGCGACCCG
CACCTGGGCCAGCGAGAGCGCAGTCGCCACGCCAGGCCCGCCAGGGGATCCTGAG**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

| | |
|--------------------|--------------|
| Restriction Sites: | Sgfl-MluI |
| ACCN: | NM_001013784 |
| Insert Size: | 900 bp |



[View online »](#)

| | |
|-------------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_001013784.1</u> , <u>NP_001013806.1</u> |
| RefSeq Size: | 4002 bp |
| RefSeq ORF: | 900 bp |
| Locus ID: | 432582 |
| Cytogenetics: | 11 B5 |