

Product datasheet for MC214945

Gm5460 (NM_001034880) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids

Tag: Tag Free

Symbol: Gm5460

Synonyms: EG432838

Mammalian Cell Selection: Neomycin

Vector: pCMV6-Entry (PS100001)

E. coli Selection: Kanamycin (25 ug/mL)

Fully Sequenced ORF: >MC214945 representing NM_001034880
Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTGCGACTCGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGCC

ATGACCACTACTCTGCCCCAGAAACCTCTACAGTTCTGACAACTTCAACAACAAACCCCTATCCCCAAAC
CAACTATTCAGGCTCCCTACTATACTCCCCCTCCTCCCGCTCCTCGTCCTCCTTCACCTGCTCC
TCCCCACTTCTCCACCTCCTCCTCGGCCACCTCCACCTCTACCTCTACCTCCACCTCCACCTGCTCC
CCAACCTTCACTTGTGCTCTGCCGCCAGCTGAGCCAAGCAAGGAAGTTTGTCTCAGCACCCTGGTACTGC
TCCTGATGCTGACCCCTGTAGTGCTGTGGTGGTCTGGTGGCTGTGCTGCAGACCGCCTGTCAAGGCCCT
GCCAGCACCTGTACCACAACCAACGCCCAAGAAAGAAAGGGCCCTCCTCCAGTGCCTCCTCCAGTG
CCTCCCCAGCTCCACCCAGCACCAGGCTCACTGGGCGCCCATTTGAATCATCTGCTGCTGTACCTGCT
ATGATGTGTGTGTGAGAAGAGGAACAGAGGGCCGTATATCAACCTTGACCTCCTGAAGGACCTCTGTTCT
TCGAGTTCCCTGCACCTCTAAGACTGGCTTCTTCCACCTGAGGAGTACCTGTGCTGGCCTCCTGTTCT
CAGTGCCACCATCCTCCACACAGTTGTCCAGGCCATTTCTTTGCTCCCTTCCCTGCCAGGTGTCTC
CCAAAGCCCCCTTGTCACTGCCTCCATGCAGCGGACCTCTCAGACACCACAGAGGCAATATGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI

ACCN: NM_001034880

Insert Size: 765 bp



OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001034880.2, NP_001030052.1</u>
RefSeq Size:	1903 bp
RefSeq ORF:	765 bp
Locus ID:	432838
Cytogenetics:	14