

Product datasheet for **MC215006**

Rhox2d (NM_001081669) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Rhox2d (NM_001081669) Mouse Untagged Clone
Tag: Tag Free
Symbol: Rhox2d
Synonyms: EG434760; Rhox2.4; Rhox2c
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC215006 representing NM_001081669
Red=Cloning site **Blue**=ORF **Orange**=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGAGCAACAAAGCATCAATTACAAGCTTGATGTGGGACCTGAGGAGGATGAGGAAAATGCGAATGGTG
TAAAGACTCTGATGGTGTGCTGGCTGGAGAGGGAAGAAACGAGGGAGAGAGTGGACCGGCCTGCCTGG
GTCGGGAGTCTCAGCAGCGGAAGGATACAGAGCAGGAGAATTAAGTGCAGGTGGGCTTGCTCGCCAGTA
GCCGACCTCATGGATAAGAGCAACCAAGAGGACCTTAGTGCCACTGGCTGTGCCAGGAGAAGGAGAAGC
AGCCAGAGGAGCCAGTCCCTGATTCCATGGGGATTTGGAAAATGTAAAGCCTATGTCAGGGCCGTGGTC
CACTGTAAATCCTGTGAGAGTGTGGTGCCGAATTCGCTACAGTTGGCAACAGAGCTTCAATGTGCTG
CAACTACAAGAGCTGGAGAGCATCTCCAGTGCAATCAATACATCAGCACTACAGAGGCAAAACGTCTGG
CAAAATCCATGGGTGTGAGTGAAGCCACAGTGCAGGAATGGTTTTTGAAGAGGAGAGAGAAAATACAGGAG
TTATAAGAGGCTG**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001081669
Insert Size: 576 bp
OTI Disclaimer:

Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001081669.2</u> , <u>NP_001075138.2</u>
RefSeq Size:	795 bp
RefSeq ORF:	576 bp
Locus ID:	434760
Cytogenetics:	X A3.3