

Product datasheet for **MC215009**

Rhox8 (NM_001004193) Mouse Untagged Clone

Product data:

| | |
|---------------------------|------------------------------------------------------------------------------------|
| Product Type: | Expression Plasmids |
| Product Name: | Rhox8 (NM_001004193) Mouse Untagged Clone |
| Tag: | Tag Free |
| Symbol: | Rhox8 |
| Synonyms: | Tox |
| Mammalian Cell Selection: | Neomycin |
| Vector: | pCMV6-Entry (PS100001) |
| E. coli Selection: | Kanamycin (25 ug/mL) |
| Fully Sequenced ORF: | >MC215009 representing NM_001004193 Red=Cloning site Blue=ORF Orange=Stop codon |

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGCTTTTCTTCTTAAGGCTTAACCAAGAGGTTAAGCCCAGAATTTCCAGAAGCAGGAGCACCGTCAACA
GCATCCATGCTATGGAACCTCAAGAAAGTCAACCCAGTCGTCTCTTCTGAGAGATGATGAAATTAAGAAAG
CGATGACGCAGCTGCATGGATAGTCTCACAGGAGATGAAGGAGAGAGAGGAGAACCGAGGTATACAGGAC
TGCCCTTTGTTGGGGCCATGGCTGCAGAAGGAGAAGGTAAAAATGAAGAAGAAGGAGAAGAAGGAGAAG
AAGGAGAAGAAGGAGAAGAAGGAGAAGAAGGAGAAGAAGGAGAAGAAGGAGAAGAAGGAGAAGAAGGA
AGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGGA
GAAGAAGGAGAAGAAGGAGAAGAAGGAGACGAAGGAGGAGAAGGAGGAGAAGGAGGAGAAGAAGGAGAAG
AAGAAGGAGAAGAAGAAGGAGGAGAAGGAGAAGGAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGA
AGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGAAGCAGCAGCAGCCGCCGCCAGAGAT
GAACTACTGCTGGCAGTTCAGCCCCAGTAGATGACAGAAGCCATGATGCCGGGGCCAGCAGCAATGGCG
AGGATAGGGGTCAAGGAGAGGAGCTGATCCCTGGGGGCACAAAGGGCCTTGAGGCTCTCCAAGTCCATC
AGGCCAGTACCCCGCAACCGCTACAGGTTTACCAAGTTCAGCTCCAGGAAGTGGAGCGCATTTTCGAA
CGCAATCACTATCCAGTGTCTGCGGCCGAAGGGAGCTCGCGAGATGGATAGGTGTGACTGAATCCAGAG
TGGAGAAGTGGTTAAGAGTAGGAGAGCCAAATACAGGAAATGCCTGATG**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja2279_h05.zip

Restriction Sites: SgfI-MluI



[View online >](#)

| | |
|-------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| ACCN: | NM_001004193 |
| Insert Size: | 963 bp |
| OTI Disclaimer: | <p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p> |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | NM_001004193.2 , NP_001004193.1 |
| RefSeq Size: | 1206 bp |
| RefSeq ORF: | 963 bp |
| Locus ID: | 434768 |
| Cytogenetics: | X A3.3 |