

Product datasheet for MC215027

Shf (NM_001013829) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Shf (NM_001013829) Mouse Untagged Clone
Tag: Tag Free
Symbol: Shf
Synonyms: RP23-173H17.6
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC215027 representing NM_001013829
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGGAGCCCTACGAGGCCAAAAGATGATGGCTGAGATCCGGGTTCTAAGGAGACTGCAGCTCAACCC
TACCTCTCTATGACACACCCTATGAGCCAGAAGATGAGGGGCCAGCCAGAAGGTGAGGGGACTCCCTG
GCCAGGGAGTCGCGCCTGCCAGAGGATGACGAGAGGCCCTCCAGAAGAGTATGACCAGCCCTGGGAGTGG
AAGAAGGAACGGATTTCCAAAGCCTTTGCTGCCAGTTTGAAGGATCGGAGAAGTGCCTGTACCTGGCC
GGGAGGAGAAGGGCGGCTACCTCCCAGGCTATCTGCTGGGAACCCAAAGTCAGCCAAGCCCTAGGCAT
GGAACCCAGCAGCCCTGGGGGAGTGGACAGATCCAGCACTGCCTCTGGAAAACCAAGTCTGGTACCAC
GGGGCCATCAGCCGAACAGATGCCGAGAACCTGCTCAGGCTGTGCAAAGAGGCCAGCTACCTGGTGC
ACAGTGAGACCAGCAAAATGACTTCTCTCTCCCTCAAGAGCAGTCAGGGCTTCATGCACATGAAGCT
GTCCCGGACTAAGGAACACAAGTATGTGTTGGGCCAGAACAGCCACCCTTCAGCAGCGTCCCTGAAATC
GTGCACCACTATGCCAGCCGAAGCTGCCATTAAGGGAGCGGAACACATGTCCCTGCTCAACCTGTGG
CCATCCGACTCTG**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001013829
Insert Size: 717 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001013829.2</u> , <u>NP_001013851.2</u>
RefSeq Size:	1564 bp
RefSeq ORF:	717 bp
Locus ID:	435684
UniProt ID:	<u>Q8CG80</u>
Cytogenetics:	2 E5
Gene Summary:	Adapter protein which may play a role in the regulation of apoptosis in response to PDGF. [UniProtKB/Swiss-Prot Function]