

Product datasheet for MC215055

Itlnb (NM_001007552) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Itlnb (NM_001007552) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Itlnb
Synonyms:	Itln1b; Itln2
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC215055 representing NM_001007552 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACCCAACGGCTTCCTGCTGTTTATCATGATTGCCACGAGAGTGTGCAGTGCAGCTGAAGAGAACC
TGGACACCAACAGATGGGGCAATCTTTCTTTCTCTGCCCAGAAGCTGTAAGGAAATCAAGCAGGA
GGACACAAAGGCACAAGATGGTCTCTATTTCTGCGCACGGAGAATGGTGTCTCTACCAGACCTTCTGT
GACATGACCACTGCAGGTGGTGGCTGGACCCTGGTGGCTAGTGTGCACGAGAACAACCTGCGTGGGAGGT
GCACTGTGGGTGATCGCTGGTCCAGTCAGCAAGGCAACAGAGCTGATTACCCAGAGGGGGATGGCAACTG
GGCCAAC TACAACACCTTTGGGTCTGCAGAGGGTGCCACAAGTGATGACTACAAGAACCCTGGCTACTTC
GACATCCAGGCAGAGAACCTGGGCATCTGGCATGTGCCCAACAACAGCCCCCTGCACACCTGGAGGAACA
GCTCCCTGCTGAGGTACCGCACCTTCACTGGCTTCCTGCAGCGCTTGGGCCATAATCTGTTTGGTCTCTA
CCAGAAGTATCCGGTGAATATGGAGAAGGAAAGTGTGGACTGACAATGGCCCAGCATTTCCTGTGGTC
TATGACTTTGGTGATGCTCAGAAGACAGCCTCTTACTCTCCCTCTGGCCGGAATGAATTCACCTGCAG
GATATGTTTCAGTTCAGAGTGTTCATAATGAGAGAGCAGCCAGTGCCTTGTGTGCTGGCGTGAGGGTCAC
TGGATGTAATACTGAACATCACTGCATCGGTGGAGGAGGATTCTCCAGAATTTGACCCCGAGGAGTGT
GGAGACTTTGCTGCATTTGATGCGAATGGATATGGAAC TACATTCCGGTACAGCAATAGCCGGGAGATAA
CTGAAGCAGCTGTGCTTCTGTTTATCGC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001007552



[View online »](#)

Insert Size:	942 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001007552.1</u> , <u>NP_001007553.1</u>
RefSeq Size:	1136 bp
RefSeq ORF:	942 bp
Locus ID:	493583
UniProt ID:	<u>Q80ZA0</u>
Cytogenetics:	1 H3
Gene Summary:	May play a protective role in the innate immune response to parasite infection. [UniProtKB/Swiss-Prot Function]