

Product datasheet for **MC215068**

Taar9 (NM_001010831) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Taar9 (NM_001010831) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Taar9
Synonyms:	Gm231; taR-9
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC215068 representing NM_001010831 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACAAGCGACTTCTCCCCAGAACCGCCCATGGAGCTCTGCTATGAGAATGTGAACGGATCCTGCATTA
AAAGCTCTTACGCACCCTGGCCTCGAGCCATCCTCTATGGGGTCTCGGTTTGGGAGCCCTGCTGGCAGT
GTTTGGGAACCTACTGGTCATCATTGCCATCCTCCACTTAAACAGCTGCACACGCCTACGAACTCCTG
GTGGCCTCCCTGGCCTGTGCTGACTTCTTGGTGGGGTGACAGTGATGCCCTTCAGCACGGTGCCTTCTG
TGGAGAGCTGCTGGTACTTTGGGAGAGTTACTGTAAGTCCACACGTGTTTCGACACCTCCTTCTGTTT
TGGCTCTCTGTTTCATTTATGCTGCATCTCCATTGACCGGTACATTGCAGTTACAGATCCACTGACCTAT
CCGACCAAGTTCCTGTATCGGTTTCTGGACTGTGCATTGCTCTCTCTTGGTTCTTTTCTGTACATACA
GCTTTTCCATCTTTTACACGGGAGCCAATGAGGAAGGGATTGAGGAATTAGTGGTTGCTCTGACCTGTGT
GGGAGGCTGCCAGGCTCCACTGAATCAAATGGGTTTTACTTTGTTTCTTTTGTCTTTCTGCCACT
GTGGTCATGGTGTCTATATGGTCGGATATTTTGGTGGCAAAGTATCAGGCTAGGAAGATAGAGGGTA
CAGCCAACCAAGCTCAGGCCTCCTCTGAGAGCTATAAGGAAAGAGTAGCCAAAAGAGAGGAAGGCTGC
CAAGACCTTGGGATCGCCATGGCTGCATTTCTCGTGTCTGGCTGCCATACATTATTGATGCGGTTATT
GACGCTTACATGAACCTCATAACTCCTGGTACGCTATGAGATATTAGTGTGGTGTGTTTACTACAATT
CAGCTATGAACCCCTTGATATACGCCTTTTTTTATCCTTGGTTTCGTAAGGCAATAAAACTATTGTGAG
TGGCAAAGTCTTCAGGGCTGATTCAACAACCTAATCTATTCTCTGAAGAGGCAGGTGCAGGTTAA

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Chromatograms: https://cdn.origene.com/chromatograms/ja2331_f02.zip



[View online »](#)

Restriction Sites:	Sgfl-Mlul
ACCN:	NM_001010831
Insert Size:	1047 bp
OTI Disclaimer:	<p>Due to the inherent nature of this plasmid, standard methods to replicate additional amounts of DNA in E. coli are highly likely to result in mutations and/or rearrangements. Therefore, OriGene does not guarantee the capability to replicate this plasmid DNA. Additional amounts of DNA can be purchased from OriGene with batch-specific, full-sequence verification at a reduced cost. Please contact our customer care team at custsupport@origene.com or by calling 301.340.3188 option 3 for pricing and delivery.</p> <p>The molecular sequence of this clone aligns with the gene accession number as a point of reference only. However, individual transcript sequences of the same gene can differ through naturally occurring variations (e.g. polymorphisms), each with its own valid existence. This clone is substantially in agreement with the reference, but a complete review of all prevailing variants is recommended prior to use. More info</p>
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	NM_001010831.1 , NP_001010831.1
RefSeq Size:	1047 bp
RefSeq ORF:	1047 bp
Locus ID:	503558
UniProt ID:	Q5QD04
Cytogenetics:	10 A4
Gene Summary:	Orphan olfactory receptor specific for trace amines.[UniProtKB/Swiss-Prot Function]