

Product datasheet for MC215156

Vmn1r170 (NM_001166722) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Vmn1r170 (NM_001166722) Mouse Untagged Clone
Tag: Tag Free
Symbol: Vmn1r170
Synonyms: EG546944; Gm5999
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC215156 representing NM_001166722
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTTGCTTCAGAATAAATCTGTGAAAACCACAGAGGAAGTGGCTCTTCAGGTCCTTTGCTTTGCCAGG
TTGGGGTTGGGACTGCAGTCAACATTTTTCTGATTGTTTCATAATTTATCTCCAACTTGAATGGATCTCA
ACAGAGACCCATAAAAGTCATTCTAGCCAACTTAGCTGTGGGCAATACCTTGATTCTCCTCTTTGCATTT
CCAAACAATATGACACTTTTTGTTCCAAAGGATCCCCAACTGACCTAAAATGTAACCTGGGTACTTCA
TTTGGCTGGTGGCTCGAAGTACAAACATGTGCTCCACCTGCTCTCTGAGCACTTATCAGCTTGTCACTCT
TGCTCCTGGTACCTGGGGTAGGGTATGCTCCGTGGAAGAGCCCCAAATTTGTGAGGTATATTTGTTAT
AGTTGTTGGTGTTCAGTGTCTTAAATAATGCTCATATCCAATGAAAGTCAGTGGTCCCAGAAAAACAC
ATAATGATACCAATACTAAAACAAATGGGTTTGCTCCACCCTGGTTTCAGTATAGGCATGCGCATCTT
GTCATTTGCCATGATGCAGTGTTCATCAGCATCATCTGGAGCAGTGTCTCCATGGTATTCTCCTG
AACAGACACCACCAGACTGCAGCATATCAATCCACAAATCAGAAGCTCAGAGTTCATGCTGAGAGTA
GAGCATCCACACCATCCTGATGCTGGTGGTCACATTTGTGACCTGTTATCATCTGGACTGATTTGTGAT
ATTCTGTAACATTTCTTTTTGGACTCTCAGCTCTGGTTGAGGCGTGTCAAACAATTTGGATGCAAGC
TTCCCTACCTTTCTCCATTACTCTTGATCTTTAGGGATCCTAAGGATCCTTGTCTCTGCTATTCTGCT
GCTGA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001166722
Insert Size: 915 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none"> 1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001166722.1</u> , <u>NP_001160194.1</u>
RefSeq Size:	915 bp
RefSeq ORF:	915 bp
Locus ID:	546944
Cytogenetics:	7 A3