

Product datasheet for MC215242

Tmem236 (NM_001081310) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Tag:	Tag Free
Symbol:	Tmem236
Synonyms:	2010003H20Rik; Fam23a; Gm613
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC215242 representing NM_001081310 Red=Cloning site Blue=ORF Orange=Stop codon TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC GCCGCGATCGCC ATGGCTTCTGGCAAAATAATCAAGCTGGTGGTTTTTGAGCTTCTTGAGTTGCTGCCTTCTCCATTCTA CCCTTGATGATAATGGAACAGTTTGGGACTGCCAATCAGCGAACAAGTCTGAAAGGACTCACTATTGGCT CATTGTTTCTGTCCATCGCCTATGTGGCTGTAGTGAGTCTCTTGATTGGGTTCCGGTAAAAGTTGTC CTATACAAGAAACGGCATCTTTACAAGAAATTATAGGATGGCGACCTGTTCTGGTGATGTGTGTGTGC TCACCACTCTCCCAAGCTTCAGCTTCTCCATAGCAGTCACCGAGGTTCAAAGAATATCAATGGCTCAGC TAATTCGTTGCCTGAATCCTTACCTGATCTGCCAGTATCTCTGGTTCTGTTGCCCTGATTGTGGTTGAC ATTATTGAAAACTCAGGCAATATCCTCTAAGAGGCAGTCAAAGGGTTATGAAGACAATGACATCTGTA TAACCTCGTTGCAACAGATAAAAACCGTGACAGAACAGGTAGTTCAGAGTGACGAAATCCTGCATCTGC CCAGGCGGCCAAGCCCACTGCGATGTCCAGCCCCGGAACCATGTGGCAGTGCTTGCCGGGCCCCCTGGAG CCCTCCTTTCACTAGGATCCTTAGGACGATGTCTCAAAGAGATGTTCCGGCGGAGCTGTTCTGCGAA GCTTCCTCATGTGGGCCGACACGGTGGAATGCTGCGTGTGGCTGGGCACCAGCGGTGTACAAGTCAGC CTGGCTGTACCCGGTCTACATATTCAGTTTTATCTCTCTCTGAGAATGGTCTTCACTCCTAAGAATCCT CTCCTTAATTCCTGGGATCTTGATGCAAGATTGCCCTTTGTTTTCTTAGACTCAGTTAATCCTTG CCCTGGGAACGATCAGGCCAGTACTGGCCTTTGTAAAAATGCTCCTAGTGACTATATCCTATATCTACTT CAATTATGTAACCAAACTCAGGCCCTTCTGCTTTGAAACCTCTCCATTTAA ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAATGATATCCTGGATT ACAAGGATGACGACGATAAGGTTTAA
Restriction Sites:	SgfI-MluI



ACCN:	NM_001081310
Insert Size:	1035 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001081310.2</u> , <u>NP_001074779.1</u>
RefSeq Size:	3658 bp
RefSeq ORF:	1035 bp
Locus ID:	625286
UniProt ID:	<u>A2ARJ3</u>
Cytogenetics:	2 A2