

Product datasheet for MC215287

Defb23 (NM_001037933) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Defb23 (NM_001037933) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Defb23
Synonyms:	EG629114
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC215287 representing NM_001037933 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGAAGCTCCTTTTCTCTCTTTGCCAGTCTCATGCTTCAGTCTCAGGTGAAATCAGAATTTGCGACCG
TGAAAAGATGCTTAGTGGGCTTTGGGAAATGCAAAGACAGCTGCCTCGCTGATGAAACCCAGATGCAGCA
TTGCAAAGCCAAAAGTCTGCATTGGCCCAAAGTGTCAGAAGTATGAAAAGCTACCTGCGACAAGAA
ATACCCACATACCGGATGAGGACATCGTAGAGATGATGAGAACGAAAAAGAATCCACAGAGGAAATGC
AAGGAAAACAAGCCTTAATGGCACTCTTCAAAGCCGAGCTGCCAAGAACCTTTGTCCAATACCAACTC
TGCTATGGGCCCAAACGCCTTCCCGTGAAGTCTGTACCACCAGGACCAGGAGACGCCGATGGGCGCT
ACTGCTTCCACCAAGATGCACATCAAACAAAGCAAGGATTCTGCCAATGCTGCCCTCGGCCCGACCTG
GACCACCT**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001037933
Insert Size:	501 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001037933.2, NP_001033022.1</u>
RefSeq Size:	624 bp
RefSeq ORF:	501 bp
Locus ID:	629114
Cytogenetics:	2 G3