

## Product datasheet for **MC215289**

### Sult2a3 (NM\_001101586) Mouse Untagged Clone

#### Product data:

**Product Type:** Expression Plasmids  
**Product Name:** Sult2a3 (NM\_001101586) Mouse Untagged Clone  
**Tag:** Tag Free  
**Symbol:** Sult2a3  
**Synonyms:** EG629203; Gm6955  
**Vector:** pCMV6-Entry (PS100001)  
**E. coli Selection:** Kanamycin (25 ug/mL)  
**Cell Selection:** Neomycin  
**Fully Sequenced ORF:** >MC215289 representing NM\_001101586  
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGGATCGCC**

ATGATGTCAGACTATTATTGGTTTGAAGGTATACCTTTCTGCCATAGCATATCAAAAAGAAATTTGG  
AAGATATTCGTAATAAGTTTGTGGTGAAGAAGAAGACTTGTGATTTGACTTACCCCAAGTCAGGAAC  
ACACTGGCTGATTGAGATTGTCTGCTTGATTTCAGACCAGGGGAGATCCAAAGTGATTCAAAGTGTGCC  
ATTTGGGACCGCTCACCTGGCTAGAGACTGATATAGGATATCTTGCATTAATTAATAAGGAAGGACCAC  
GACTCATTTCCTCTCATCTTCTGTCCATCTCTTCAAGTCTTTTTTTCAGTTCCAAGGCAAGGCGAT  
CTATCTCATCAGAAATCCAGAGATATTCTGTGTCTGGTTACTTTTTCTATGGTAACACAAACCTGTG  
AAGAATCCAGGATCACTCAGAATATTTTGAATGGTTCCTCAAAGGAAATGTCATATATGGATCATGGT  
TTGAGCATGTCCGTGGCTGGCTGTCCATGAGAGAATGGGATAACTTCTGGTACTGTACTATGAAGACAT  
GAAAAAGGACACAAAGGGAACATAAAGAAGATCTGTGACTTCTAGGGAAAAAATTAGAACCAGATGAG  
CTGGATCTCGTCTCAAATACAGCTCTTCCAAGCCATGAAAGAAAACAACATGTCCAATTATAGTCTCG  
TTAGCGAAGATATGATTACTAATGGCTTGAAGCTCATGAGAAAAGGCACAACCTGGGGACTGGAAGAATCA  
TTTCACAGTAGCCAAGCTGAAGCCTTCAATAAAGTTTTCCAGGAGAAAATGGCTGGATCCCCCAAGGG  
ATGTTCCCATGGGA**TAA**

**ACGCGT**ACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

**Restriction Sites:** SgfI-MluI  
**ACCN:** NM\_001101586  
**Insert Size:** 858 bp



[View online »](#)

<b>OTI Disclaimer:</b>	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u><a href="#">NM_001101586.2</a></u> , <u><a href="#">NP_001095056.2</a></u>
<b>RefSeq Size:</b>	858 bp
<b>RefSeq ORF:</b>	858 bp
<b>Locus ID:</b>	629203
<b>Cytogenetics:</b>	7 A1
<b>Gene Summary:</b>	This is one of seven sulfotransferase family 2A genes in a chromosome 7 A1 cluster. [provided by RefSeq, May 2010]