

## Product datasheet for **MC215350**

### Cypt9 (NM\_001039942) Mouse Untagged Clone

#### Product data:

|                           |                                                                                    |
|---------------------------|------------------------------------------------------------------------------------|
| Product Type:             | Expression Plasmids                                                                |
| Product Name:             | Cypt9 (NM_001039942) Mouse Untagged Clone                                          |
| Tag:                      | Tag Free                                                                           |
| Symbol:                   | Cypt9                                                                              |
| Mammalian Cell Selection: | Neomycin                                                                           |
| Vector:                   | pCMV6-Entry (PS100001)                                                             |
| E. coli Selection:        | Kanamycin (25 ug/mL)                                                               |
| Fully Sequenced ORF:      | >MC215350 representing NM_001039942<br>Red=Cloning site Blue=ORF Orange=Stop codon |

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
 GCC**CGCATCGCC**

ATGACTCGGGTGGCCAAGAAAGTCCACTGGTCTAGAGCAGCAACAGCAGCAGCCGCCGCTGCCAAAGCCA  
 AAAGATCGAAGCTTAAGAAGAAGCGAAACCCTAGATCTAAGCTTCCCAAGAGATCTCGTCACTCCTTAAT  
 TCGTTCTCATTCTCGTTCCTGTTGTTGCTGTCTTCGTTCTCGTAGCAGAACTACATTAAGATTACGAAG  
 AAGGGGAACAATCACTGCGAAGGAGAATTAGGAGACGAATCAAGCGTGCTACTGAACCTTCGATTGATGC  
 GAAGTCAAATGGAGCAGAGCCAACGGGAGCTGATCGAACCAGGAACTGAACTTTGGAGCCAAATGA  
 GATTACAGTGGCATTGTTCTCTCATGAGAATGCCAATGTGTCTGATCCAGAGGAAGTTCCACCATGCCTT  
 GACAGTGACCCATTTCCGAATGGAGACTTGCCAGTTCC**TAG**

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
 ACAAGGATGACGACGATAAGGTTTAA

|                    |                                                                                                                                                                                                                                                                                                                                                |
|--------------------|------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Restriction Sites: | SgfI-MluI                                                                                                                                                                                                                                                                                                                                      |
| ACCN:              | NM_001039942                                                                                                                                                                                                                                                                                                                                   |
| Insert Size:       | 462 bp                                                                                                                                                                                                                                                                                                                                         |
| OTI Disclaimer:    | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |


[View online »](#)

|                        |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
|------------------------|-------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| Components:            | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).                                                                                                                                                                                                                                                                                                                              |
| Reconstitution Method: | <ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol> |
| Note:                  | Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.                                                                                                                                                                                                                                                                                                                                                                                          |
| RefSeq:                | <u>NM_001039942.2, NP_001035031.1</u>                                                                                                                                                                                                                                                                                                                                                                                                                                                                             |
| RefSeq Size:           | 533 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| RefSeq ORF:            | 462 bp                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Locus ID:              | 664724                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                            |
| Gene Summary:          | DISCONTINUED: This Gene record was discontinued due to lack of supporting evidence.                                                                                                                                                                                                                                                                                                                                                                                                                               |