

Product datasheet for MC215365

Mthfd2l (NM_026788) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Mthfd2l (NM_026788) Mouse Untagged Clone
Tag: Tag Free
Symbol: Mthfd2l
Synonyms: 1110019K23Rik; C630010D07Rik
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC215365 representing NM_026788
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGGATCGCC**

ATGGCGACGCGGGCCCGGGCTTCTCGCTGCTCCGGGGCCGTCTGGGCCGAGGGCCTGTGCGAGCCCCGG
GCGTCGCCGAGCGCGCTGGCGCGGCTTCGGGAGCAGCGCCGGAGACACGAAGCCGTTATCATCTCAGG
AACTAACATGGCCAAACAGATCCAGAAGGAGATACAGCAAGGAGTGAATCATGGATCGCTCTCGGGAAC
AGGAGGCCCATCTGAGTATCATTTTAGTGGGAGATAACCCAGCAAGCCACACCTATGTCAGGAATAAGA
TAAAAGCGGCTTCTGCTGTAGGTATCTGCAGTGAGCTCATATAAAGCCCAAGAATGTTTCTCAGGAGGA
ACTGTTGGACATCACTGATCAGCTGAACATGGATCCAAGAGTCAGCGGTATATTAGTTACAGTTCCACTC
CCAGACCATGTAGATGAGAGAACAATA TGCAATGGAATTGCTCCAGAAAAGGACGTTGATGGATTCACA
TTATCAATATTGGACGACTGTGCCTTGACCAGCATTCACTGATCCCTGCCACTGCCAGTGGCCTTTGGGA
AATCATCAAAGAGCAGGAATCGAAACATTTGGGAAAAATGTGGTCGTGGCGGGAAGATCCAAGAACGTG
GGCATGCCGATTGCCATGCTGCTGCACACGGACGGAGAGCAGCAGCGTCTGGAGGTGATGCAACGGTGA
CGATAGCGCACAGATACACCCAGAGAACTGAAGGCCACACTCAACTAGCAGACATCATCATAGT
GGCTGCGGGTATCCAGGTTGATTACAGCTGACATGGTCAGAGAAGGTGCGGCAGTAATTGATGTGGGT
ATCAACTATATAAGATCCTGTGACGGAAAAGACCAAGTTAGTTGGAGATGTGGACTTTGAAGCTGTTA
AGAAGAAAGCCAGCTTCACTCCAGTCCCAGGAGGTGTGGACCCATGACCGTGGCCATGCTTCTTAA
GAACACACTTTTAGCTGCAAAGAATATCGCCTAC**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_026788



[View online »](#)

Insert Size:	1017 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_026788.1</u> , <u>NP_081064.1</u>
RefSeq Size:	2256 bp
RefSeq ORF:	1017 bp
Locus ID:	665563
UniProt ID:	<u>D3YZG8</u>
Cytogenetics:	5 E1