

Product datasheet for **MC215442**

Gm11567 (NM_001101613) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gm11567 (NM_001101613) Mouse Untagged Clone
Tag: Tag Free
Symbol: Gm11567
Synonyms: OTTMUSG00000002206
Mammalian Cell Selection: Neomycin
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Fully Sequenced ORF: >MC215442 representing NM_001101613
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACCAACTCCTGCTGCTCCCCTTGCTGCCAGCCACCTGCTGCAGGACCACCTGCTGCAGGACCACCT
GCTGGAGACCAAGCTGTGTGACCAGCTGTTGCCAGCCCTGCTGCCAGCCAGCTGCTGTGGCTCCAGCTG
CTGCCAGCCTTGCTGTCAAACCACCTGCTGCAGGACCTGCTTCCAGCCATGCTGTGTGAGCAGCTGCTGC
CGCACACCCTGCTGCCAGCCCTGCTGCTGTGTGCCAGCTGCTGCCAGCCCTGCTGCCAGCCAGCTGTT
GTCAGTCCAGCTGTTGCCAGCCAGCTGCTGTCAGTCCAGCTGCTGCCAGCCTAGCTGTTGTCAGTCCAG
CTGCTGTCAGCCAGGTGCTGTATCTCCAGCTGCTGCCAGCCCTGCTGCAGGCCTAGCTGCTGTCAGTCC
AGCTGCTGCAAGCCCTGCTGCCAGCCCTTCTGCCTCAACCTGTGCTGCCAGCCAGCCTGCTCGGACCTG
TGACCTGCACCAGGACTTGCTACCAGCCAACATGTGTCTGTGTGCTGCTGCCAGCCTGCTGCCAGGCTGTGG
GTCTAACTGCTCCCAATCTTGCTGT**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: Sgfl-Mlul
ACCN: NM_001101613
Insert Size: 588 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001101613.1, NP_001095083.1</u>
RefSeq Size:	1043 bp
RefSeq ORF:	588 bp
Locus ID:	670533
Cytogenetics:	11