

Product datasheet for **MC215455**

Zglp1 (NM_001103168) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Zglp1 (NM_001103168) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Zglp1
Synonyms:	Glp1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC215455 representing NM_001103168 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTAAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGAGGCGGCACAGGCTGGAGACTTGACGCGCCGGCAGGAATTCTGGCGCCACCTTGCTGGACACCG
AGTCGCTCCGGAAGAGTCGACCCCCGCGCTGGAGCCTGGAGCCCTGAGATGCCTCACTCCCAATATCAG
GTCTCTGTGGCTACATGCCAGGACTCGGTGAGCACTGCATTGCCTTTCCTCAAGAGAAAGAAAAGGGA
CTGCCGGGTTCCCTTCGCCGGCCACCCAGGTCTGGGGTCGTGCTGGGAGCTGATGGTCATAGGAATGT
CAGACCACCTGAGTATGGCCAGGAATCCCAGGGGCACACAGTGCCTAACCTGGAATTTCCAGTGCTAC
TTCTCCTGCCTCACTCCAGAGAAGACCACGGAAGCAGCTAAACCCGCGAATGGGTATTGAGAAGGTGGAC
CCCAGGTTCAAGGGGTAACCTCTGGAGTTCCAGATCCAGCCGGATTCTAGCCTCCAGATTGTACCCACCT
ACAGCTTGCTGGAAGAAGTTGTTCAAAAAGCTTCTGCAAGCCCTTCCAAGGCTCTGGCCAGCCAGG
AAGCAGTGAGGCCCTTGGGCCAGGCGGTGTGCTTCTGTAGAACCAGAGGACTCCACTATGGAGAGAC
GCTGAAGATGGGACCCCTCTGCAATGCCTGTGGCATCAGGTACAAGAAATATGGCACCCGTTGTTCCA
GTTGCTGGCTGGTCCCAAGGAAAAGCATCCAGCCCAAAGGCTCTGTGGCAGGTGTGGGATGTCCAGGA
CCCCACCTCAGCCCAACCCAGGAAGT**TAA**

ACGCGTACGCGGCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001103168
Insert Size:	801 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001103168.1</u> , <u>NP_001096638.1</u>
RefSeq Size:	1010 bp
RefSeq ORF:	801 bp
Locus ID:	100009600
UniProt ID:	<u>Q1WG82</u>
Cytogenetics:	9 A3
Gene Summary:	Transcriptional repressor that plays a central role in somatic cells of the gonad and is required for germ cell development. Able to repress GATA transcription factor function. [UniProtKB/Swiss-Prot Function]