

Product datasheet for **MC215457**

Tmsb15b2 (NM_001080967) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Tmsb15b2 (NM_001080967) Mouse Untagged Clone
Tag: Tag Free
Symbol: Tmsb15b2
Synonyms: Tb15c
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC215457 representing NM_001080967
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAGCGATAAACCAGACTTGT**CAGAAGTTGAAACATTTGACAAAGCAAAGTTGAAGAAA**CTA**ACTG**
AAGTAAAAATA**CTCTTCCGTCGAAGGAACTATCCAGCAAGAGAAAGAACA**CAAT**GAAAGAACA****TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001080967
Insert Size: 138 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).



[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

RefSeq: [NM_001080967.4](#), [NP_001074436.1](#)

RefSeq Size: 691 bp

RefSeq ORF: 138 bp

Locus ID: 100034363

Cytogenetics: X

Gene Summary: This gene is part of tandem duplication; the other copy of this gene is about 1.5 kb away. Readthrough transcripts that include two exons transcribed from each gene exist. [provided by RefSeq, Jul 2008]