

Product datasheet for MC215522

Gm2663 (NM_001102660) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Gm2663 (NM_001102660) Mouse Untagged Clone
Tag: Tag Free
Symbol: Gm2663
Synonyms: 100040208
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC215522 representing NM_001102660
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGAAGATCATCTTCTTCTTCACTTTTCTGGGAGCTGCTGTTGCTCTCCAGCTAACAGTGATGACAAGA
TTGTTGGAGGCTACACATGTCCGAAGCATTTCGTTTCCTTACCAGGTGTCTCTGAATGATGGCATTAGCCA
TCAGTGTGGTGGCTCCCTTATCAATGATCAGTGGTACTGTCTGCTGCTCACTGCTACAAAAGAAGACTC
CAGGTTTCGCCTTGGTGAACACAATATTGATGTCCTTGAGGGTGGTGGAGCAATTTATTGATGCAGAAAAGA
TCATTCGACACCCAGACTATAACAAGGACACTGTGGATAATGACATCATGCTGATTAATGAAGTCACC
TGCCATCCTTAACTCCCAAGTATCCACAGTCTCTCTGCCAGATCCTGTGCATCTACAAATGCTCAGTGC
CTTGTTTCTGGCTGGGAAATACTGTGAGCATTGGTGGAAATACCCAGCACTCCTTCAATGTCTGGAAG
CTCCTGTCTCTCGGCCAGTCTTGCAAAAAATCCTACCCAGGCCAGATTACCAGCAATATGTTCTGCCT
GGGATTCCTGGAGGGTGGAAAGGACTCCTGTGATGGTATTCTGGTGGACCTGTTGTCTGCAATGGAGAG
ATCCAGGGTATTGTCTCCTGGGGTTCAGTCTGTGCAATGAGAGGGAAGCCTGGTGTTCACACCAAGTCT
GCAACTACCTAAGCTGGATTTCAGGAGACCATGGCGAACAAC**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001102660
Insert Size: 744 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001102660.1</u> , <u>NP_001096130.1</u>
RefSeq Size:	870 bp
RefSeq ORF:	744 bp
Locus ID:	100040208
Cytogenetics:	6