

## Product datasheet for MC215525

### Gm11563 (NM\_001126320) Mouse Untagged Clone

#### Product data:

Product Type:	Expression Plasmids
Product Name:	Gm11563 (NM_001126320) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Gm11563
Synonyms:	2310045F10Rik; OTTMUSG00000002200
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC215525 representing NM_001126320 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC  
GCC**CGATCGCC**

ATGGTCAGCTCCTGTTGTGGCTCTGTCTGCTCTGAGGAGGGCTGTGGCCAAGGCTGCTGCCAGCCAGCT  
GCTGCCAGACCCTGCTGTAGGACCACCTGCTGCCGCCCCAGCTGTTGTGTGTCCAGCTGCTGCAGGCC  
TCAGTGTGCCAGTCTGTGTGTGCCAGCCTACCTGCTGCCGCCCCAGCTGTTGCCGCCCCAGCTGCTGT  
GTGTCCAGCTGCTGCAGACCTAGCTGCTGCCAGTCTGTGTGCTGCCAGCCACCTGCTGTCGCCCCAGCT  
GCTGCCGCCCCAGCTGCTGCATCTCCAGCTGCTGTCAACCCTCCTGTGGTAGCTCCAGCTGTTGTGGTTC  
CAGCTGCTGCCGCCCCCTGTTGCCGCCCCCTGCTGCCGCCCCCTGCTGCTGTTTGTGAGACCAGTCTGTGGTCAG  
GTCTGCTGCCAAACCACTTGCTACCGCCCACTTGCTGTCATCTCCACCTGCCCCGCCCCATGTGCTGTG  
CCACCCCTGCTGCTGA

**ACGCGT**ACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT  
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001126320
Insert Size:	507 bp
OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

<b>Components:</b>	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
<b>Reconstitution Method:</b>	<ol style="list-style-type: none"><li>1. Centrifuge at 5,000xg for 5min.</li><li>2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.</li><li>3. Close the tube and incubate for 10 minutes at room temperature.</li><li>4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.</li><li>5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.</li></ol>
<b>RefSeq:</b>	<u>NM_001126320.2, NP_001119792.1</u>
<b>RefSeq Size:</b>	1018 bp
<b>RefSeq ORF:</b>	507 bp
<b>Locus ID:</b>	100040248
<b>Cytogenetics:</b>	11 D