

Product datasheet for MC215538

Gm16430 (NM_001166601) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gm16430 (NM_001166601) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Gm16430
Synonyms:	EG666120
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC215538 representing NM_001166601 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**C

ATGGCTCTTAAGAACTGTGGCGATACCAAAGGACGGTTACTTATTACTTCTTGACTATGATGATGAAG
ACGATGACATAAATTTTTGGAGGAGGCTCATTCTGAAGAAAACGTAAGTTTCTCAGAGGAATGGCAGCG
TTTTGCAAGTTCTGTAGAGACACCAATAGAAAACAGGAATTTGTTAAGTGGTGAACAGCAAGATGGGAAT
GCTTCAAATTTGGATCTTATGGAAGAACAGAATCCAGTAACTCATGATGATGAGAATGAAATTCCTGAAG
AAATCGTTGGAGATACACGGGAGATGATCAACAACAAGTCGTGTGAGCAGTATAAACTACATTTTCAGAA
GTTTGATATGGATGTCCAGAATTTCAATGAACAACAAGAAAAATCAGTGGGTTTGATGAACCTGGAGACC
AACAACCTCGATATGCTTTTGTAGTAGATGGTGAATTGAGAAAA**TGA**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001166601
Insert Size:	468 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).


[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	1. Centrifuge at 5,000xg for 5min. 2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA. 3. Close the tube and incubate for 10 minutes at room temperature. 4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom. 5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
Note:	Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.
RefSeq:	<u>NM_001166601.1, NP_001160073.1</u>
RefSeq Size:	700 bp
RefSeq ORF:	468 bp
Locus ID:	100040894
Cytogenetics:	X A5