

Product datasheet for MC215579

Vmn1r79 (NM_001166835) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Vmn1r79 (NM_001166835) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Vmn1r79
Synonyms:	Gm9807
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Cell Selection:	Neomycin
Fully Sequenced ORF:	>MC215579 representing NM_001166835 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTGAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCCGCGATCGC

ATGAAATGGATCCTGGAAATTTGCAATAGAACTTTCTCTTCTGCCAGATTACAGTGGGAATGTTTG
GCAATTCCTCAATCTATTTTATTATATTATTTGATATTTACAGAAAAGCATTTAACGCCCAAAGACCT
GATTATAGAGCACTTGACTTTCCCAACTACTTGACTGTCATCTTAAGAGGGTTCCACAGACAATGTCA
TATTTTGGGTTTAAAGATTTCTAGATGACATGGGATGTAATTAATAGTGATATTTCCCGGATAACAA
GGGGGGTGTCTTTATGCTATGTGCCTATTGAGTTGTTTCCAAGCAATCACCATCAGTCCCAATAACTC
CAGGTGGATAAAGCTTAAATACAGAGCAACCAAGTACATTGGTCCCTCCTGCTCAGTCAGTTGGCTTG
CACCTGCTTCTAAACATCTTGACTCCACTAAGAGAATCAGGCACCAATTACAAGAAAAATGCAACTAACA
GGTTGAGTTATGGATACTGCTCATTGTTTGCTTCTGACAGTGTAGCAACTACATTGTATATATCTTAGT
GTGCTTCTCTGATGTCTTCTGTCTGGGTCTCATGGCCTGCTCAAGTATCTCCATGGTTAGTATTCTCTAC
AGACATAAGAGACAAGTCAAGCATATCCACAGTGTCAACACTTTCTGAAAGTCTCACCTGAAGACAGAG
CCACCCAACTATCCTCATCTGCTGTGCATTTGTCTCTTACTCATTCTCTCCATTGTGGCTAT
CATTAGAACCTACTTAACTATCCAGTGATATGGGGAGTAAATATATTTACATCCCTAGAAATATTCTTT
CCCATGTTTTGTCCCTTTGTTCTCATCACAATATGAGATCTAGTTCCAATCTGTTTTTACCATGCTTAG
GTAAGAGATAA

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001166835
Insert Size:	921 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001166835.1</u> , <u>NP_001160307.1</u>
RefSeq Size:	921 bp
RefSeq ORF:	921 bp
Locus ID:	100042437
Cytogenetics:	7 A1