

Product datasheet for MC215591

Gm10267 (NM_001165950) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gm10267 (NM_001165950) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Gm10267
Synonyms:	ENSMUSG00000069385; VtIgl1
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC215591 representing NM_001165950 Red=Cloning site Blue=ORF Orange=Stop codon

TTTGTATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGC**

ATGTTTCATCTTCTCGAGAGTTCAAGTTATATATATTACCTTATCCTTTCTTTGTTTTCAGAAAGTCACT
TTATAAGAAGAGCAATTTACAAGCACATAGTAATGTGTAGAGGGTTAGCAGTAAAGAATTTGTACCAG
AGAGTATTTTCCAGTCTGTGCAACCAATGGCCGTACTTATTTCAACAAATGCATTTCTGTTTAGCCTAC
AGAGAAAATGATGGTTCATTTATTATGTCACATCTCGGCAAATGC**TAA**

ACGCGTACGCGGCCGCTCGAGCAGAAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	SgfI-MluI
ACCN:	NM_001165950
Insert Size:	258 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).

Components: The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).


[View online »](#)

Reconstitution Method:

1. Centrifuge at 5,000xg for 5min.
2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.
3. Close the tube and incubate for 10 minutes at room temperature.
4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.
5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.

Note: Plasmids are not sterile. For experiments where strict sterility is required, filtration with 0.22um filter is required.

RefSeq: NM_001165950.1, NP_001159422.1

RefSeq Size: 258 bp

RefSeq ORF: 258 bp

Locus ID: 100042855

Cytogenetics: 18 B3