

Product datasheet for MC215629

Gm4307 (NM_001166641) Mouse Untagged Clone

Product data:

Product Type:	Expression Plasmids
Product Name:	Gm4307 (NM_001166641) Mouse Untagged Clone
Tag:	Tag Free
Symbol:	Gm4307
Synonyms:	100043239
Mammalian Cell Selection:	Neomycin
Vector:	pCMV6-Entry (PS100001)
E. coli Selection:	Kanamycin (25 ug/mL)
Fully Sequenced ORF:	>MC215629 representing NM_001166641 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGACTACCATGGCTGACAAAATGGACATGTCATTGGAAGACATCATTAAAGCTGAACAAGATGCAGCAAG
GACGTCGTGACAGGCCTGATCCCGGTGAAGCGAGGCACTGGCCTAAGAGGTATCGGCCAGCTTTCAC
CCATGGTGAAGGAACAGGCTGGCACCATACTGCAGACCAAAACAACCTCCAGACAAGTGGCAGCATGAC
CTATTCATTGGAGGCTTCAGGGTCAAACCACGTGGATACTGGTGGGAAGCTGTTCTGTCAAATCTGC
ACTTCGGAGTGTGAGATGCTGATATTCAGCTACTCTTTGCTGAATTTGGAACGTTGAAAAATCTGCTGT
GCACTATGATCGCTGTGGACGAAGTTTAGGGACAGCACATGTGCACTTTGAAAGGAAAGCAGATGCCCTG
AAGGCTATGAGAGAGTACAATGGCGCCCTTTGGATGGCCGCCCTATGAACATCCACCTTGTCACTTCAC
AGATTGATAGACAAGGAAGACCTGCACTAAACAGTGATAAGGGTGGCATGACAAGAAACCCTGGCTCTGG
AGTATTAAGTGGTGGAGGCACCAAGAGATGGCACTTGGAGGCAGCCAGGGAAGCGGGAGAGGCACCAGC
AGGAACTCAAACCTGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGCAGC
AACAGCAGAAGCAGCAGCAACAGCAGCAGAAGCAGGGACAGAATCACCAGCACCAGCAGCAGCAGAAGGA
GCAGCAGCAGCAACAGCAGCAGAAGGAAGTGCAGCAGCTATCTGTAGAGGAGCTGGACGCCACGCTGGAT
TATTATCAGGATATGATGGACACCAGCTAA

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites:	Sgfl-MluI
ACCN:	NM_001166641
Insert Size:	870 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001166641.1</u> , <u>NP_001160113.1</u>
RefSeq Size:	910 bp
RefSeq ORF:	870 bp
Locus ID:	100043239
Cytogenetics:	10