

Product datasheet for **MC215635**

Vmn1r100 (NM_001166844) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Vmn1r100 (NM_001166844) Mouse Untagged Clone
Tag: Tag Free
Symbol: Vmn1r100
Synonyms: ENSMUSG00000074331; Gm10671
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC215635 representing NM_001166844
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**GCGATCGCC**

ATGTCTGCTCATGGTAAATCCGTGAAAACCACTGAGGAAGTGGCTCTTCAGATCCTTTGCTTTGCCAGT
TTGGGATTGGGACCGTGGCCAATGTCTTCTGTTGTCCATAATTTCTCTCCAGTCTTGACTGGTCTAA
ACAGAGGCCAGACAGGTGATTTAAGCCACATGGCTGTGGCCAATGCCCTTGACTCTATTCCTCACTATA
TTTCCAAACAACATGAGTGCTTTTGTCTCCAAAACCTCCTCCAACCTGAACTCAAATGTAAATTAGAATTCT
TCAGTCACATGGTGGCAAGAAGCACAAATTTGTGTTCCACCTGTGTCTGAGTGTCCATCAGTTTGTCCAC
TCTTGTTCCTGTTAATAGAGGTAAGGTAAGTAACTTATACTCAGAGCAAGTGTCCAAAATTTGCGAATTAT
TCTTGCTACAGTTGTTGGTTTTTCAGTGTCTTAAGTAACATCCACATTCCAATTAAGGTCAGTGTCCAC
AGATAATAGACAATAACACTGACTCTAAAAGCAACTTGTCTGTTCCACTTCTGGATTATTGTAGGCAT
TGTCTTCTGCAGTTTTCCCATGATGCCACATTCATGAGCATCATGGTCTGGACCAGTGTCTCCATGGTA
CTTCTCCTCCATAGACATCGCCAGCGAATGCAGCACATCCTCACTCCAATCAGGATGCCAGAGGCCAAG
CTGAGACCAGAGCAACCCATACTATCCTGATGTTGGTGGTCACATTTGTTAGCTTTTATCTTCTAAATTG
TATTTGTATCATCTTTCATGCCTTTTCTATACATTCTCGTCTCTTCATAAGGCTTGTGAGGTTCTG
GCTGCAGTCTTCCCAGTATCTGCCCTTACTGTTGATCTTCAGAGATCCTAAAGATCCTTGTCTGTGC
TCTTCAAATGTTGA

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001166844
Insert Size: 924 bp



[View online »](#)

| | |
|-------------------------------|---|
| OTI Disclaimer: | Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP). |
| Components: | The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water). |
| Reconstitution Method: | <ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C. |
| RefSeq: | <u>NM_001166844.1</u> , <u>NP_001160316.1</u> |
| RefSeq Size: | 924 bp |
| RefSeq ORF: | 924 bp |
| Locus ID: | 100043536 |
| Cytogenetics: | 7 A3 |