

Product datasheet for MC215642

Vmn1r137 (NM_001166849) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: Vmn1r137 (NM_001166849) Mouse Untagged Clone
Tag: Tag Free
Symbol: Vmn1r137
Synonyms: 100043614; Gm4551
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC215642 representing NM_001166849
 Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
 GCC**CGGATCGCC**

ATGTCTGCTCATGGTAAATCCGTGAAAACCACTGAGGAAGTGGCTCTCCAGATCCTTTGCTTTGCCAGT
 TTGGGGTTGGAAGTGTGGCCAATGTCTTTCTGTTGTCCATAATTTCTCTCCAGTCTTGACTGGTTCTAA
 ACAGAGGCCAGACAGGTGATTTAAGCCACATGGCTGTGGCCAATGCCTTGACTCTATTCCTCACTATA
 TTTCCAAACAACATGAGTGCTTTTGTCTCCAAAACCTCTCCAAAATGAACTCAAATGTAAATTAGATTCT
 TCAGTCACATGGTGGCAAGAAGCAAAAATTTGTGTTCCACCTGTGTCTGAGTATCCATCAGTTTGTCCAC
 TCTTGTTCCTGTTAATAGAGGTAAGGTAAGTAACTTATACTCAGAGCAAGTGTCCAAAATTTGTGGAATTAT
 TCCTGCTACAGTTGTTGGTTTTACAGTGTCTTAAGTAACATCCACATTCCAATTAAGGTCAGTGGTCCAC
 AGATAACAGACAATAACACTGACTCTAAAAGCAACTTGTCTGTTCCACTTCTGGATTCAATGTAGGCAT
 GGTCTTCTGCAGTTTTCCCATGATGCCACATTCATGAGCATCATGGTCTGGACCAGTGTCTCCATGGTA
 CTTCTCCTCCATAGACATCGCCAGCGAATGCAGCACATCCTCACTCCCAATCAGGATGCCAGAGGCCAAG
 CTGAGACCAGAGCAACCCATACTATCCTGATGCTGGTGGTGCACATTTGTTAGCTTTTATCTTCTAAATTG
 TATTTGTATCATCTTTCATGCCTTTTCTATACATTCTCGTCTCTCATAAGGCTTGTGAGGTTCTG
 GCTGCAGTCTTCCCAGTATCTGCCCTTACTGTTGATCTTCAGAGATCCTAAAGATCCTTGTCTGTGC
 TCTTCAAATGTTGA

ACGCGTACGCGGCCGCTCGAGCAGAAAACCTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
 ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001166849
Insert Size: 924 bp



[View online »](#)

OTI Disclaimer:	Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).
Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001166849.1</u> , <u>NP_001160321.1</u>
RefSeq Size:	924 bp
RefSeq ORF:	924 bp
Locus ID:	100043614
Cytogenetics:	7