

Product datasheet for **MC215654**

1700047117Rik2 (NM_001100116) Mouse Untagged Clone

Product data:

Product Type: Expression Plasmids
Product Name: 1700047117Rik2 (NM_001100116) Mouse Untagged Clone
Tag: Tag Free
Symbol: 1700047117Rik2
Synonyms: MGC61221
Vector: pCMV6-Entry (PS100001)
E. coli Selection: Kanamycin (25 ug/mL)
Cell Selection: Neomycin
Fully Sequenced ORF: >MC215654 representing NM_001100116
Red=Cloning site Blue=ORF Orange=Stop codon

TTTTGTAATACGACTCACTATAGGGCGGCCGGGAATTCGTCGACTGGATCCGGTACCGAGGAGATCTGCC
GCC**CGATCGCC**

ATGGAGCGTGAGTCGGGATGCGCGCGGGGAGAGACCGAGGCAGCGGCCGCCACGGCCTTCCGAGACG
CGACACGGCAGATAAGTAATGAAAGAGGATTTGAGAATGTAGAACTGGGCGTCATGGGAAGAAGAAGAA
AGTCCCGAGGAGGTCATCCACTTTGTGCTGAGGAGACAATGGAGGAGTACAGCACCGACGAGGACGAG
GTGGACGGCCTGGACAAGAAGGACGTGCTGCCACTGTGACCCGACAAAATTACTGGGGCCCTACC
TATGGTTTTACATGCTGCGGGCCGCCACATCAACCCTCTCAGTGTGTGACTTTCTTGAGAGAAGATTGC
GTCAGTTTTGGGCATCAGCACCCAAAATACAGTACGCCATTGATGAGTATTACCGGATGAAGAAGGAG
GAGGAGGAGGAGGAAGAAGAAAAAGGATGTCAGAAGAAGCAGAAAGACAATACCAACAAAAAGCTGC
AGGCCGATTCCATTGTACAGACAGACCAACCAGAAACAGTGTGCTCCAGCTTTGTAATATTAATTTTGA
AATGGAGGAAGACTGTGAAGCAATTAAGGAAAACAAACAGCGTCCAGTCTCTGTCCACCA**TAG**

ACGCGTACGCGGCCGCTCGAGCAGAACTCATCTCAGAAGAGGATCTGGCAGCAAATGATATCCTGGATT
ACAAGGATGACGACGATAAGGTTTAA

Restriction Sites: SgfI-MluI
ACCN: NM_001100116
Insert Size: 624 bp

OTI Disclaimer: Our molecular clone sequence data has been matched to the reference identifier above as a point of reference. Note that the complete sequence of our molecular clones may differ from the sequence published for this corresponding reference, e.g., by representing an alternative RNA splicing form or single nucleotide polymorphism (SNP).



[View online »](#)

Components:	The ORF clone is ion-exchange column purified and shipped in a 2D barcoded Matrix tube containing 10ug of transfection-ready, dried plasmid DNA (reconstitute with 100 ul of water).
Reconstitution Method:	<ol style="list-style-type: none">1. Centrifuge at 5,000xg for 5min.2. Carefully open the tube and add 100ul of sterile water to dissolve the DNA.3. Close the tube and incubate for 10 minutes at room temperature.4. Briefly vortex the tube and then do a quick spin (less than 5000xg) to concentrate the liquid at the bottom.5. Store the suspended plasmid at -20°C. The DNA is stable for at least one year from date of shipping when stored at -20°C.
RefSeq:	<u>NM_001100116.1</u> , <u>NP_001093586.1</u>
RefSeq Size:	3863 bp
RefSeq ORF:	624 bp
Locus ID:	100101807
UniProt ID:	<u>Q8BR63</u>
Cytogenetics:	12 C1